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LiquiBand® FIX8™ Hernia Mesh Adhesive Fixation Device

Single use Instrument

Read all directions, precautions and warnings prior to use.

These instructions for use provide directions for the proper use of the Advanced Medical Solutions (Plymouth) Ltd. LiquiBand® FIX8™ Hernia Mesh Adhesive Fixation (HMF) device. They are not intended to be a comprehensive instruction manual for the performance of a laparoscopic surgical procedure. Prior to the use of the device, surgeons should be fully proficient in minimally invasive/laparoscopic and open hernia repair surgical techniques.

Device Description

LiquiBand® FIX8™ HMF device is designed for the application of n-butyl-2-cyanoacrylate adhesive to an implanted hernia repair mesh, in order to fix the mesh to the underlying tissue. The device consists of:

- n-butyl-2-cyanoacrylate adhesive monomer, in liquid form, supplied in a thin-walled, sealed glass vial; and,
- a surgically invasive, laparoscopic delivery instrument comprising a 5mm diameter cannula, with a handle at the proximal end incorporating a loading chamber, filter, piston chamber and trigger. The distal tip of the device is open to allow the adhesive to be dispensed from it.

The device is designed for introduction and use through a 5 mm diameter laparoscopic port sleeve. A larger diameter sleeve will require the use of a converter.

Both the cyanoacrylate adhesive in the glass vial and the surgically invasive delivery device are supplied sterile, for single use only.

Principle of Operation

The glass vial containing the liquid cyanoacrylate adhesive monomer comes preloaded into, and is subsequently broken in, a loading chamber in the handle of the delivery instrument. The adhesive is drawn through a filter into a piston chamber in the delivery instrument handle. After priming, each press of the trigger on the handle of the delivery instrument dispenses 12.5 mg of adhesive from the distal tip of the cannula. A single device contains sufficient adhesive for at least 33 individual applications of the adhesive. The gauge on the side of the delivery instrument gives an indication of the amount of adhesive delivered and the approximate amount remaining.

Mechanism of Action

When applied to the proximal surface of the hernia mesh, the liquid adhesive monomer permeates through the perforations in the mesh to the surface of the underlying tissue, where it polymerises on contact with moisture on the tissue surface. This process of chemical polymerisation fixes the mesh to the surface of the tissue at the site of adhesive contact, maintaining the mesh in position while it is incorporated into the abdominal wall through the normal process of tissue fibrosis.

The adhesive completes its polymerisation reaction within approximately 10 seconds. Once completely set, the adhesive no longer possesses adhesive properties such that tissues or surgical instruments may be placed in contact with it without risk of unwanted adhesion.

Indications for use

The LiquiBand® FIX8™ HMF device is intended for use in the laparoscopic surgical repair of inguinal hernia, achieved through the fixation of polypropylene or polypropylene/polyester combination hernia mesh to the abdominal wall.

Contraindications

- The device is not intended for use when prosthetic material fixation is contraindicated.
- Do not use on patients with a hypersensitivity to cyanoacrylate adhesives or formaldehyde.
- Do not use for the fixation of meshes constructed with polytetrafluoroethylene (PTFE) or absorbable materials.
- Do not use device for closure or fixation of cerebral tissues, blood vessels or peripheral nerves.

Warnings

- Only physicians having adequate training and familiarity with endoscopic techniques should use this device. A thorough understanding of the operating principles, risks versus benefits, and hazards involved in utilizing an endoscopic approach is necessary to avoid possible hazards to the user and/or patient.
- This device is provided sterile and is intended for use in a single patient. Do not reuse, reprocess, clean, disinfect or resterilize this device as this may compromise the sterility and performance of the device.
- Ensure the device is properly activated and primed before use according to the instructions in the Directions for Use section.
- This device does not contain any user-serviceable parts. Do not attempt to repair or dismantle the device. If at any point the device appears to be damaged or not functioning correctly, discard and replace it with another device.
- Do not dilute or mix the adhesive with other substances.
- Accidental bonding of unwanted tissue may occur due to misapplication of adhesive. Separation of tissues after accidental bonding should only be performed if deemed necessary.

Precautions

- The fixation method for any prosthetic device should be determined on the basis of accepted surgical techniques, procedural requirements, and the instruction for use of the prosthetic material.
- Ensure the mesh is held in contact with the underlying tissue during each application of the adhesive and for approximately 10 seconds afterwards, to allow for polymerization.
- The viscosity of the adhesive is only slightly greater than that of water, so adhesive should be applied very carefully to prevent its spread to unwanted areas.
- The adhesive should always be applied in minimal amounts, i.e. avoid multiple applications of adhesive in any given location. A second application of adhesive can be applied over the first only after full polymerization.
- The application of an excessive amount of adhesive in a single location prolongs setting time and can prevent adherence. After polymerization, any excess adhesive may lead to detachment of the adhesive film and/or give rise to the formation of small fragments of polymerised adhesive.


- The adhesive will readily adhere to most substrates. Care should be taken to avoid unwanted contact with the adhesive during polymerization. Polymerized adhesive can be removed from metal instruments with acetone.
- There is only limited medium to long term clinical data available for the safe use of this device for the closure of the peritoneal sac.
- There is only limited clinical data to support the safe use of this device for the fixation of polyester hernia repair meshes. Bench tests have demonstrated suitable compatibility and performance.

Potential Adverse Reactions

As with the majority of implanted devices, adverse reactions associated with the use of this device may include transient local irritation at the implant site and a transitory inflammatory foreign body response.

Directions for Use

Priming the device

The device must be properly activated and primed to ensure that it operates correctly.

- Remove the device from its packaging. Hold the device with the handle vertical and the distal tip of the cannula tilted downwards.

Lift the blue lid firmly upward until a cracking sound of the glass ampoule breaking inside the lid is heard.

- Wait for 5 seconds then, whilst the tip pointing downwards and handle held vertically slowly pull out red tab no. 2, observing that the adhesive is being drawn into the transparent barrel. This should take approximately 20 seconds. At the end of piston travel gently disengage the red tab and dispose of.

Warning: A fast transfer stage could result in excessive air being drawn into the dispensing chamber.

- Gently close the blue lid until it is completely aligned with the body and clicks into place.
- Remove the red pull tab no. 4 and dispose of.
- Rotate the blue knob clockwise approximately 320° until blue dispenser trigger alongside the handle is released. Care should be taken to avoid impeding the release of the trigger. The distal tip should be directed into the packing tray well no. 6 before the blue knob is rotated.
- With the cannula tip placed into packing tray well no. 6, prime the device by actuating the trigger, approximately three times, until a droplet of adhesive freely emerges.
- Wipe any excess adhesive from the tip of the cannula onto the inside of the packing tray or on a sterile wipe. The device is now ready for surgical use.

Using the Device

- Adhesive should be applied in discrete locations distributed across the mesh to fix the mesh in place.
- Upon complete fixation of the mesh, surgical graspers should be used to test for the adhesion of the mesh and for any areas left unadhered.

Disposal

After use, dispose of the device according to local procedures and guidelines.

Storage

This device should always be stored in its original packaging. Store only at temperatures between 5°C (41°F) and 25°C (77°F). Do not use the device after the expiration date shown on the blister pack.

Sterility

This device is packaged for single patient use and is sterilized by irradiation and ethylene oxide gas. Do not re-sterilize. Do not use if package is open or damaged.

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