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**URGENT**

**VOLUNTARY ADVISORY NOTICE 015  
CONCERNING THE USE OF  
GAMBRO PRISMAFLEX™ CONTINUOUS RENAL REPLACEMENT SYSTEM**

Gambro Renal Products is voluntarily issuing an Advisory Notice for the Gambro Prismaflex™ Continuous Renal Replacement System – Model No. 6023014700.

The Advisory Notice gives instructions towards important issues in relation to the use of the Prismaflex™ equipment:

1. **Prismaflex™ Service Mode**

Gambro Renal Products is releasing together with this Advisory Notice a Calibration Instruction re-emphasising the importance of correct procedures when performing Scale and Syringe calibration.

Gambro Renal Products wants to emphasize the risk of unauthorised access to the Service Mode.

Any access given to the Prismaflex™ Service Mode by other personnel might lead to patient injury or death in case the Service Mode procedures are incorrectly carried out (e.g. Scale Calibrations).

In order to diagnose and solve potential technical problems, access to the Prismaflex™ Service mode should only be provided to trained and authorised personnel. The Prismaflex™ Service Mode is for that reason password protected.

**IMPORTANT NOTE:**

The enclosed attachments below must be made available to the trained and authorised personnel responsible for calibration of the Prismaflex machines.

The new attached Spare Parts Instruction SPI should be exchanged with the old version that you currently have.

When you place an order for the Calibration Weight Kit as a spare part you will automatically receive the new updated Calibration Instructions.

Gambro Renal Products apologizes for any inconvenience this may cause.

If you have any questions, please call your local Gambro representative.

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## **CALIBRATION INSTRUCTIONS**

### **Prismaflex™**

#### **Instructions for Scale Calibration**

When calibrating the Prismaflex™ scales in Service mode:

The potential risk of a fluid imbalance of the patient exists if the incorrect calibration weights are used during the scale calibration procedure.

As an immediate action an improved Spare Part Instruction (SPI) is released together with this Calibration Instruction that clearly explains how to assemble the Prismaflex™ calibration weight kit (G5000101) when used for calibration of the scales.

The following warning is added in the SPI:

**“Warning:**

When calibrating the scales, for the second reference point (5,200 Kg), you must use the correct weights (Part A + Part B = 5,200 Kg). Failure to use the correct calibration weights for the second reference point can cause serious injury or death to the patient”.

#### **Instructions for Syringe Pump Calibration**

When calibrating the Prismaflex™ syringe pump in Service mode:

A situation of potential false syringe alarms, e.g. “Syringe clamp” or “Empty” of the syringe pump can occur if the incorrect reference weights are used during the syringe pump calibration procedure.

As an immediate action an improved Spare Parts Instruction (SPI) is released together with this Calibration Instruction that clearly explains how to assemble the Prismaflex™ calibration weight kit (G5000101) when used for calibration of the syringe pump.

The following caution is added in the SPI:

**“Caution:**

When calibrating the syringe pump, for the second reference point (7,000 Kg), you must use the correct weights (Part A + Part B+ Part C+ Part D = 7,000 Kg). Failure to use the correct calibration weights for the second reference point will set the over-pressure alarm too low for the syringe”.

**Note:** In the Prismaflex™ Operator’s Manual Chapter “Before you get started Precautions, section 5” it is clearly stated that the accuracy of the Prismaflex™ Control Unit depends on accurate scale and pressure calibration. Ensure that scales and pressure sensors are accurately calibrated. Calibrations **must** be performed by a trained and qualified person.

Jan Spjuth  
Intensive Care Customer Support  
Gambro Lundia AB  
Monitor Division

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**HOSPAL S.A.S.**