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12" @ 100%

 **Dade Behring Inc.**  
Newark, DE 19714, U.S.A.

**EC REP** **Dade Behring Limited**  
Regus House, Atterbury  
Milton Keynes MK10 9RG  
United Kingdom

2005-03 E PN 755423.001

**Dimension® PBNP Insert Sheet**  
PN 755423.001 Rev E  
Issued 2004-03

6.0" @ 100 %

# DADE BEHRING

**Dimension®** clinical chemistry system  
**Heterogeneous Immunoassay Module**

## Flex® reagent cartridge

## PBNP

### N-terminal Pro-Brain Natriuretic Peptide Flex® reagent cartridge

**Intended Use:** The PBNP assay used on the Dimension® clinical chemistry system with the heterogeneous immunoassay module is an *in vitro* diagnostic assay for the quantitative determination of N-terminal pro-brain natriuretic peptide (NT-proBNP) in human plasma. In individuals suspected of having congestive heart failure (CHF), measurements of NT-proBNP are used as an aid in the diagnosis and assessment of severity. The test is further indicated for the risk stratification of patients with acute coronary syndrome and heart failure.

**Summary:** Left ventricular dysfunction can occur as part of coronary heart disease, arterial hypertension, valvular disease, and primary myocardial disease. If the left ventricular dysfunction remains untreated and is progressive, the potential for mortality is high, e.g., due to sudden cardiac death.

Chronic cardiac insufficiency is a clinical syndrome caused by impairment of the cardiac pumping function. Based on the symptoms, the severity of cardiac insufficiency is classified in stages I-IV by the New York Heart Association (NYHA). Clinical tests and imaging procedures are used to diagnose left ventricular dysfunction.<sup>1</sup>

The significance of natriuretic peptides in the control of cardiovascular system function has been established. Initial studies reveal that natriuretic peptides can be used for diagnosing clinical problems associated with left ventricular dysfunction.<sup>2</sup> The following natriuretic peptides have been reported in the literature: atrial natriuretic peptide (ANP), brain natriuretic peptide (BNP), and C-type natriuretic peptide (CNP).<sup>3,4</sup>

ANP and BNP have natriuretic and diuretic properties. As antagonists of the renin-angiotensin-aldosterone system, they influence the electrolyte and fluid balance in an organism.<sup>5,6</sup> In subjects with left ventricular dysfunction, serum and plasma concentrations of BNP increase, as do the concentrations of the biologically inactive prohormone, proBNP. ProBNP consists of 108 amino acids. It is secreted mainly by the left ventricle of the heart and, in this process, is cleaved into physiologically active BNP (77-108), and the N-terminal fragment NT-proBNP (1-76).<sup>4</sup>

Studies indicate that NT-proBNP can be used in diagnostic and prognostic applications.<sup>7-13</sup> The concentration of NT-proBNP in plasma indicates the prognosis for the left ventricular dysfunction. It is also useful in assigning symptoms to cardiac or non-cardiac causes.

Determination of NT-proBNP helps to identify subjects with left ventricular dysfunction. Changes in NT-proBNP concentration can be used to evaluate the success of treatment in patients with left ventricular dysfunction. There are indications that NT-proBNP, due to its functions, is suitable for use in assessing vascular remodeling, and therefore contributes to the establishment of individualized rehabilitation procedures.<sup>10,11</sup>



Fisher et al found that congestive heart failure patients with elevated NT-proBNP values had a one year mortality rate of 53% compared to 11% in patients with lower values<sup>12</sup>. In the Gusto IV study which involved more than 6800 patients it was shown that NT-proBNP was the strongest independent predictor of one year mortality in patients with acute coronary syndrome.<sup>13</sup>

The European Society of Cardiology Task Force for the Diagnosis and Treatment of Chronic Heart Failure recommends in their guidelines that natriuretic peptides including NT-proBNP "may be most useful clinically as a rule out test due to consistent and very high negative predictive values."

NT-proBNP levels are increased in patients with unstable angina and following myocardial infarction<sup>14,15</sup>. Although not diagnostic for these conditions, several studies indicate NT-proBNP measurements provide prognostic information for the short- and long-term risk stratification of patients with unstable angina or myocardial infarction.<sup>15,16,17</sup>

**Principles of Procedure:** The PBNP method is a one step enzyme immunoassay based on the "sandwich" principle. Sample is incubated with chromium dioxide particles coated with polyclonal antibodies which recognize epitopes located in the N-terminal part of proBNP, and a conjugate reagent [alkaline phosphatase (ALP)] labeled polyclonal antibody specific for a second independent epitope on NT-proBNP, to form a particle/NT-proBNP/conjugate sandwich. Unbound conjugate is removed by magnetic separation and washing. After separation and washing, the particle/NT-proBNP/conjugate sandwich is transferred to the cuvette where the sandwich bound ALP triggers an amplification cascade.\* ALP dephosphorylates synthetic flavin adenine dinucleotide phosphate (FADP) to produce FAD. FAD binds to apo D-amino acid oxidase and converts it to active holo D-amino acid oxidase. Each molecule of holo D-amino acid oxidase produces multiple molecules of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). H<sub>2</sub>O<sub>2</sub> in the presence of horseradish peroxidase (HRP), converts 3,5-dichloro-2-hydroxybenzenesulfonic acid (DCHBS) and 4-aminoantipyrene (4-AAP) to a colored product that absorbs at 510 nm. The color change measured is directly proportional to the concentration of NT-proBNP present in the patient sample.

\*Technology licensed from London Biotechnology, Ltd., London, U.K.

### Reagents

Wells <sup>a</sup>	Form	Ingredient	Concentration <sup>b</sup>	Source
1,2	Tablets <sup>c</sup>	FADP	0.056 mg/mL	
		D-proline	68.8 mg/mL	
3,4	Tablets <sup>c</sup>	Apo-D-amino acid oxidase	0.23 mg/mL	Porcine Kidney
		Horseradish peroxidase	0.19 mg/mL	Horseradish
		DCHBS	3.9 mg/mL	
		4-AAP	0.94 mg/mL	
5	Empty			
6	Tablets <sup>c</sup>	Antibody-CrO <sub>2</sub>	1.3 mg/mL <sup>d</sup>	Sheep
7	Liquid	Buffer and stabilizers		
8	Liquid	Antibody-Alkaline phosphatase <sup>d</sup>		Sheep
		Buffer and stabilizers		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value in hydrated cartridge.

c. Tablets contain excipients, buffers, and stabilizers.

d. Antibody titer and conjugate activity vary from lot to lot.

**Reagent Preparation:** Mixing and diluting are automatically performed by the instrument.

**Precautions:** Irritant. Contains 2-chloroacetamide. May cause sensitization by skin contact. Avoid contact with skin. Wear suitable gloves. Contains sodium azide (<0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

**Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.**

For *in vitro* diagnostic use

**Storage Instructions:** Store at 2 – 8 °C.

**Expiration:** Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days. Once wells 1, 2, 3 and 4 have been entered by the instrument, they are stable for 2.5 days. Once wells 6, 7 and 8 have been entered by the instrument, they are stable for 5 days.

\***Specimen Collection:** Plasma samples (lithium heparin and sodium heparin) can be collected using normal procedures.<sup>18</sup>

Comparison testing results on the Dimension® system for lithium heparin plasma samples versus sodium heparin plasma samples is provided as follows:

Comparative Specimen	Slope	Intercept	Correlation Coefficient	N
Na Heparin	0.95	0.9	0.998	55

Model equation by Passing-Bablok linear regression is: results for Dimension® system = [slope x (Li heparin results)] + intercept.

Serum samples should not be used with the PBNP method.

Samples are stable for three days when stored at 20 – 25 °C or 2 – 8 °C or for 12 months when frozen at -20 °C or below.<sup>19</sup> Repetitive freezing and thawing of specimens should be avoided.

Ensure that patient samples, calibrators and controls are equilibrated at ambient temperature (22 - 28 °C) before testing. Samples containing precipitates must be centrifuged before performing the assay.

Blood collected in the presence of EDTA, citrate or oxalate and frozen plasma samples with insufficient anticoagulants potentially can cause clumping of the chrome particles and should not be used.

Each laboratory should determine the acceptability of its own blood collection tubes. Variations in these products may exist between manufacturers and, at times, from lot-to-lot.

### Potential Interfering Substances

The following substances have no significant effect (less than 10 %) when added to approximately 125 pg/mL NT-proBNP at the concentrations indicated:

Compound	Concentration
Acetaminophen	20 mg/dL [1323 µmol/L] <sup>e</sup>
Allopurinol	2.5 mg/dL [184 µmol/L]
Amiodarone	20 µg/mL [31 µmol/L]
Amikacin	15 mg/dL [256 µmol/L]
Amidopine Besylate	4 µg/mL [7.1 µmol/L]
Ampicillin	5 mg/dL [143 µmol/L]
Ascorbic Acid	3 mg/dL [170.3 µmol/L]

Atenolol	1.0 mg/dL	[37.6 µmol/L]
Bilirubin, conjugated	60 mg/dL	[1026 µmol/L]
Bilirubin, unconjugated	20 mg/dL	[342 µmol/L]
Caffeine	10 mg/dL	[515 µmol/L]
Captopril	5.0 mg/dL	[230 µmol/L]
Carbamazepine	12 mg/dL	[508 µmol/L]
Chloramphenicol	25 mg/dL	[774 µmol/L]
Chlordiazepoxide	2 mg/dL	[67 µmol/L]
Chlorpromazine	5 mg/dL	[157 µmol/L]
Cholesterol	500 mg/dL	[12.9 mmol/L]
Cimetidine	10 mg/dL	[396 µmol/L]
Creatinine	30 mg/dL	[2652 µmol/L]
Cyclosporine A	4000 ng/mL	[3.3 µmol/L]
Dextran75	2500 mg/dL	[333 µmol/L]
Diazepam	2 mg/dL	[70 µmol/L]
Diclofenac	60 µg/mL	[188.6 µmol/L]
Digitoxin	60 ng/mL	[78.6 nmol/L]
Digoxin	5 ng/mL	[6.4 nmol/L]
Diltiazem	120 µg/mL	[266.1 µmol/L]
Dipyridamole	30 µg/mL	[59.4 µmol/L]
Erythromycin	20 mg/dL	[273 µmol/L]
Ethanol	350 mg/dL	[76 mmol/L]
Ethosuximide	30 mg/dL	[2125 µmol/L]
Furosemide	2 mg/dL	[61 µmol/L]
Gentamicin	12 mg/dL	[251 µmol/L]
Hemoglobin (monomer)	1000 mg/dL	[0.62 mmol/L]
Heparin, sodium	8 IU/mL	[8000 IU/L]
Hydralazine	20 µg/mL	[101.7 µmol/L]
Hydrochlorothiazide	20 µg/mL	[67.2 µmol/L]
Ibuprofen	40 mg/dL	[1939 µmol/L]
Indomethacin	16 µg/mL	[36.9 µmol/L]
Isosorbide dinitrate	6.0 mg/dL	[254 µmol/L]
Lidocaine	6 mg/dL	[256 µmol/L]
Lipemia (triglyceride)	3000 mg/dL	[33.9 mmol/L]
Lisinopril	16 µg/mL	[36.2 µmol/L]
Lithium chloride	3.5 mg/dL	[5.07 mmol/L]
Lovastatin	16 µg/mL	[39.6 µmol/L]
L-thyroxine	60 µg/dL	[0.77 µmol/L]
Methyldopa	2.5 mg/dL	[118 µmol/L]
Nicotine	2 mg/dL	[123 µmol/L]
Nifedipine	6.0 mg/dL	[173 µmol/L]
Nitrofurantoin	40 µg/mL	[168 µmol/L]
Nitroglycerin	0.16 µg/mL	[0.577 µmol/L]
Oxazepam	12 µg/mL	[41.8 µmol/L]
Oxytetracycline	100 µg/mL	[217.2 µmol/L]
Penicillin G	25 U/mL	[25,000 U/L]
Pentobarbital	10 mg/dL	[442 µmol/L]
Phenobarbital	15 mg/dL	[646 µmol/L]
Phenytoin	10 mg/dL	[396 µmol/L]
Primidone	10 mg/dL	[458 µmol/L]
Probenecid	200 µg/mL	[700.9 µmol/L]
Propoxyphene	0.4 mg/dL	[12 µmol/L]
Propranolol	0.15 mg/dL	[5.1 µmol/L]
Protein, human albumin	6 g/dL	[60 g/L]
Protein, human IgG	6 g/dL	[60 g/L]
Protein, total high	12 g/dL	[120 g/L]
Protein, total low	3-4 g/dL	[30 - 40 g/L]
Quinidine	20 µg/mL	[61.6 µmol/L]
Rheumatoid Factors	500 IU/mL	[500,000 IU/L]
Salicylic acid	50 mg/dL	[3.6 mmol/L]
Simvastatin	32 µg/mL	[76.5 µmol/L]
Sulfamethoxazole	320 µg/mL	[1263.3 µmol/L]
Theophylline	25 mg/dL	[1388 µmol/L]
Trimethoprim	64 µg/mL	[220.7 µmol/L]
Urea	500 mg/dL	[83.3 mmol/L]
Uric acid	20 mg/dL	[1.2 mmol/L]
Valproic acid	50 mg/dL	[3467 µmol/L]
Verapamil	16 mg/dL	[0.33 µmol/L]
Warfarin	4 µg/mL	[13.0 µmol/L]

e. Systeme International d'Unites (SI Units) are in brackets

#### Procedure

##### Materials provided

PBNP Flex® reagent cartridge, Cat. No. RF423

##### Materials needed but not provided

Dimension® PBNP Calibrator, Cat. No. RC423

Reaction Vessels, Cat. No. RXV1A

Chemistry Wash, Cat. No. RD701

Probe Cleaner, Cat. No. RD702

Sample Probe Cleaner, Cat. No. RD703

Sample Diluent, Cat. 791092901

Quality Control products such as:

MAS® and Cardiolmmune® proBNP Liquid Assayed Cardiac Marker Control

BioRad Liquichek™ Cardiac Markers Control LT

Audit™ MicroFD™ BNP Control

MAS® and Cardiolmmune® are registered trademarks of Medical Analysis Systems, Inc. Camarillo, CA

Liquichek™ is a trademark of BioRad, Bio-Rad Laboratories, Diagnostics Group, Irvine, California 92618

Audit™ and MicroFD™ are trademarks of Audit™ MicroControls, Inc. Las Vegas, NV 89102

##### Test Steps

Sampling, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system with heterogenous immunoassay module. For details of this processing, refer to your Dimension® system manual.

##### Test Conditions

##### Reaction Vessel

Sample Volume:	50 µL
Antibody-CrO <sub>2</sub> :	25 µL
Antibody-ALP:	40 µL
Temperature	37 °C ± 1.0 °C
Incubation period	6.4 minutes

<b>Cuvette</b>	<b>Reaction</b>	<b>Blanking</b>
Transfer Volume	40 µL	0 µL
FADP Reagent Volume	30 µL	30 µL
APO Reagent Volume	30 µL	30 µL
Diluent Volume	350 µL	390 µL
Temperature	37 °C ± 0.1 °C	
Reaction Period	5.4 minutes	
Wavelength	510 and 700 nm	
Type of measurement	Bichromatic Rate	
Units	pg/mL	

##### Calibration

The general calibration procedure is described in your Dimension® system manual. The following information should be considered when calibrating the PBNP method:

**PBNP requires lot specific calibration scalars. See PBNP Flex® carton label. These scalars must be entered in the calibration screen prior to calibration. Failure to enter these scalars will cause inaccurate results.**

Assay Range	10 – 30,000 pg/mL
Reference Material:	Dimension® PBNP Calibrator (Cat. No. RC423)
Calibration Range:	0 - 30,000 pg/mL
Calibration Levels	0, 300, 1500, 12,000, 30,000 pg/mL
Calibration Scheme:	Five levels. Level 1 is tested in quadruplicate. Levels 2, 3 and 5 are tested in triplicate. Level 4 is tested in duplicate. Assure that enough sample is dispensed in each cup to allow for the number of determinations for each level during calibration.
Calibration Frequency:	Every new reagent lot. or Every 30 days for any one lot.
Starting Coefficients:	C <sub>0</sub> -2826.0 C <sub>1</sub> 16014.0 C <sub>2</sub> -2.2 C <sub>3</sub> 18293.0 C <sub>4</sub> 0.5
Calibration Scalars	See carton label. Flex® lot specific scalars must be entered in the calibration screen prior to calibration.
A, B, C and D:	

##### Quality Control

At least once each day of use, analyze two levels of quality control material with known NT-proBNP concentration. Test results are invalid if they fall outside of the laboratory's acceptable limits. Follow the procedure outlined in your Dimension® system manual.

A system malfunction may exist if the following 5-test precision is observed:

<b>Concentration</b>	<b>S.D.</b>
150 pg/mL	>15 pg/mL
450 pg/mL	>45 pg/mL

**Results:** The Dimension® system automatically calculates concentrations of NT-proBNP in pg/mL using the calculation scheme illustrated in your Dimension® system manual. Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

##### Limitations of Procedure

Results < 10.0 pg/mL should be reported as "less than 10.0 pg/mL" instead of the numerical value.

Manual dilutions: > 30,000 pg/mL  
Results in excess of 30,000 pg/mL should be repeated after diluting the sample with Dimension® PBNP calibrator Level 1 or Dimension® Sample Diluent to obtain results within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Autodilution is not recommended for this method.

One-step immunometric assays are susceptible to a high dose "hook effect", where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.<sup>20</sup> Such samples must be diluted and reassayed prior to reporting the results (see dilution information listed above). The PBNP method shows no hook effect up to at least 300,000 pg/mL.

Concentrations of natriuretic peptides may be elevated in patients with acute myocardial infarction, patients that are candidates for renal dialysis, and patients that have undergone renal dialysis.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies.<sup>21</sup>

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® system manual.

**Expected Values:** NT-proBNP concentrations in the Reference Group are shown in the following tables. Each laboratory should establish its own reference interval for the PBNP method as determined on the Dimension® system. The recommended decision thresholds from these distributions are:

Patients < 75 years:	125 pg/mL
Patients ≥75 years:	450 pg/mL

1st fold

2nd fold

12" @ 100%

### Reference Study Group

NT-proBNP concentration was determined in 308 individuals without congestive heart failure (163 women and 145 men). This population included apparently healthy individuals and individuals with diabetes, hypertension, and pulmonary disease. The statistics for NT-proBNP concentrations in the reference study group are shown in the following table.

All	<55 yrs	55 - 64 yrs	65 - 74 yrs	≥ 75 yrs
Mean	32.1	39.1	83.1	353.6
SD	38.2	46.2	54.6	775.8
Median	17.9	18.8	59.5	125.6
95 <sup>th</sup> Percentile	81.9	163.1	202.6	1372.8
% < 125 pg/mL	97%	93%	88%	-
% < 450 pg/mL	-	-	-	83%
N	163	15	8	122
Males	<55 yrs	55 - 64 yrs	65 - 74 yrs	≥ 75 yrs
Mean	29.7	37.9	76.0	414.3
SD	42.5	41.6	-	889.0
Median	13.4	13.4	76.0	111.8
95 <sup>th</sup> Percentile	134.4	114.3	76.0	1475.7
% < 125 pg/mL	95%	100%	100%	-
% < 450 pg/mL	-	-	-	77%
N	76	6	1	62
Females	<55 yrs	55 - 64 yrs	65 - 74 yrs	≥ 75 yrs
Mean	34.2	39.9	84.1	291.0
SD	34.0	51.6	58.9	639.7
Median	23.7	18.8	59.5	131.3
95 <sup>th</sup> Percentile	75.8	163.1	202.6	1080.7
% < 125 pg/mL	99%	89%	86%	-
% < 450 pg/mL	-	-	-	88%
N	87	9	7	60

### Disease Study Group

Blood samples were obtained from 227 patients diagnosed with congestive heart failure (CHF). The population included 69 women and 158 men. The descriptive statistics and New York Heart Association (NYHA) functional classes are provided below. Each laboratory should establish a reference range that represents the patient population that is to be evaluated. In addition, laboratories should be aware of their institution's current practice for the evaluation of CHF.

CHF Population - All	<55 yrs	55 - 64 yrs	65 - 74 yrs	≥ 75 yrs
Mean	2400.1	3120.7	4167.4	4579.3
SD	3855.9	8809.1	8550.7	8721.2
Median	660.5	533.3	1506.7	2513.4
95 <sup>th</sup> Percentile	10961.2	16817.2	10889.3	11398.6
% > 125 pg/mL	79%	73%	90%	-
% > 450 pg/mL	-	-	-	91%
N	48	49	61	69

CHF Population - Males	<55 yrs	55 - 64 yrs	65 - 74 yrs	≥ 75 yrs
Mean	2748.2	3475.1	4427.8	5773.0
SD	4179.3	10332.0	9669.3	10238.5
Median	1148.2	821.7	1273.4	3086.6
95 <sup>th</sup> Percentile	12127.1	16817.2	10889.3	14118.8
% > 125 pg/mL	77%	84%	89%	-
% > 450 pg/mL	-	-	-	92%
N	34	32	45	47

CHF Population - Females	<55 yrs	55 - 64 yrs	65 - 74 yrs	≥ 75 yrs
Mean	1554.7	2453.7	3435.2	2028.9
SD	2886.2	5024.5	4175.0	2581.3
Median	257.8	492.9	1828.4	1077.1
95 <sup>th</sup> Percentile	10961.2	20855.8	13287.1	8891.2
% > 125 pg/mL	86%	53%	94%	-
% > 450 pg/mL	-	-	-	91%
N	14	17	16	22

CHF Population - All NYHA Functional Class	NYHA Functional Class				
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	3693.0	1492.9	2124.5	5423.2	7134.5
SD	7911.8	2260.2	2946.1	10692.5	11900.1
Median	1422.0	659.2	1077.1	2636.5	2754.5
5 <sup>th</sup> Percentile	59.3	47.9	29.4	117.4	76.4
95 <sup>th</sup> Percentile	11398.6	5748.8	8344.4	13287.1	20855.8
% > Cutoff	85%	76%	83%	90%	90%
Minimum	17.6	21.7	17.6	60.3	50.6
Maximum	70025.3	10613.6	16817.2	70025.3	63515.3
N	227	53	72	71	31


CHF Population - Males NYHA Functional Class	NYHA Functional Class				
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	4273.6	1706.9	2582.3	6273.3	9300.3
SD	9109.3	2242.9	3238.0	12099.8	15619.8
Median	1710.7	938.6	1273.4	3269.9	2717.1
5 <sup>th</sup> Percentile	60.3	57.8	33.1	122.1	76.4
95 <sup>th</sup> Percentile	13600.1	5748.8	8344.4	15686.3	63515.3
% > Cutoff	86%	81%	83%	93%	88%
Minimum	21.7	21.7	27.3	60.3	76.4
Maximum	70025.3	10613.6	16817.2	70025.3	63515.3
N	158	41	48	53	16

CHF Population - Females NYHA Functional Class	NYHA Functional Class				
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	2363.5	761.6	1208.8	2920.1	4824.3
SD	3725.4	1301.7	2013.8	3749.6	5565.2
Median	933.5	232.6	550.2	1933.5	2754.5
5 <sup>th</sup> Percentile	50.6	47.9	24.4	103.2	50.6
95 <sup>th</sup> Percentile	10375.0	4723.5	2524.6	13287.1	20855.8
% > Cutoff	81%	58%	83%	83%	93%
Minimum	17.6	47.9	17.6	103.2	50.6
Maximum	20855.8	4723.5	10029.4	13287.1	20855.8
N	69	12	24	18	15


These results show that there is a relationship between the severity of the clinical signs and symptoms of CHF and the median NT-proBNP concentration, demonstrating that the Dimension® PBNP Method can be used as an aid in the diagnosis of all degrees of CHF severity including asymptomatic patients.

### Interpretation of Results

The Receiver Operator Characteristics (ROC) Curve presents the clinical sensitivity and specificity at various cutoffs for the 227 patients diagnosed with CHF and 308 subjects without CHF. The ROC curve for the Dimension PBNP assay is shown below. The area under ROC curve (AUC) for the Dimension PBNP assay is 0.921 with a 95% confidence interval of 0.898 to 0.943.

A box and whiskers plot of the clinical study population is presented below. Recommended clinical thresholds are 125 pg/mL for patients younger than 75 years and 450 pg/mL for patients 75 years and older. Three disease group samples with values above the assay range are not displayed in the plot.

### Sensitivity and Specificity vs. Age and Gender

The tables below show the clinical sensitivity and specificity of the Dimension® PBNP assay using a cutoff of 125 pg/mL for patients younger than 75 years and 450 pg/mL for patients 75 years or older.

Males	<75 yrs	≥ 75 yrs
Sensitivity	84% (93/111)	92% (43/47)
95% Confidence Interval	77 - 91	84 - 99
Specificity	95% (79/83)	77% (48/62)
95% Confidence Interval	91 - 100	67 - 88

Females	<75 yrs	≥ 75 yrs
Sensitivity	77% (36/47)	91% (20/22)
95% Confidence Interval	64 - 89	79 - 100
Specificity	97% (100/103)	88% (53/60)
95% Confidence Interval	94 - 100	80 - 96

### Age -matched ROC analysis

An age-matched ROC analysis of the clinical data was performed via the weighted method described in Kondratovich (2002)<sup>22</sup>. The resulting AUC is 0.925 with a 95% confidence interval of 0.903 to 0.946.

### Specific Performance Characteristics<sup>f</sup>

Material	Reproducibility <sup>g</sup>				
	Mean (pg/mL)	Within-Run Precision		Total Precision	
		SD (pg/mL)	%CV	SD (pg/mL)	%CV
Human Plasma					
Pool 1	159.0	3.4	2.2	9.1	5.7
Internal QC Pools					
Pool 1	449.5	8.1	1.8	16.6	3.7
Pool 2	956.7	15.3	1.6	34.9	3.6
Audit™ MicroFD™ BNP Control					
Level 1	175.5	2.0	1.1	6.8	3.8
Level 2	3733.8	71.8	1.9	115.9	3.1

f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (Refer to your Dimension® system manual).

g. Reproducibility testing was done in accordance with the NCCLS Approved Guideline for User Evaluation of Precision Performance of Clinical Chemistry Devices EP5-A, 1999. Specimens at each level were analyzed in duplicate once per day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

### Correlation

A method comparison study was performed to demonstrate equivalency of the Dimension® PBNP and Elecsys® proBNP assays. Samples for the study were obtained from patients known to have CHF (n=224), patients with signs and symptoms suggestive of CHF (n=17), and randomly selected patients without knowledge of medical history or diagnosis (n=111).

### Regression Statistics<sup>h,i,j</sup>

Comparative Method	Slope	Intercept pg/mL	Correlation Coefficient	n
Elecsys® proBNP	0.90 (95% CI: 0.88 to 0.92)	-15.4 (95% CI: -25.0 to -8.9)	0.985	352

h. Model equation for regression statistics is: results for Dimension® system = [Slope x (comparative method results)] + intercept.

i. The range of PBNP values in the correlation study was: 16.1 - 29,893.1 pg/mL.

j. Elecsys® Systems is a registered trademark of Roche Diagnostics GmbH, D-68298 Mannheim. US Distributor: Roche Diagnostics Corporation, Indianapolis, IN 46250-0457 USA.

### Functional Sensitivity

Functional sensitivity of Dimension® PBNP is ≤ 30.0 pg/mL. Functional sensitivity is the lowest NT-proBNP concentration that can be measured with a total precision coefficient of variation of 20%.

### Analytical Specificity

The following substances have no significant cross-reactivity (less than 1 %) at the concentrations indicated when added to samples containing 0 and approximately 125 pg/mL NT-proBNP:

Substance	Concentration
ANP <sub>28</sub>	3.1 µg/mL
NT-proANP <sub>1-30</sub> (preproANP <sub>26-55</sub> )	3.5 µg/mL
NT-proANP <sub>31-67</sub> (preproANP <sub>56-92</sub> )	1.0 ng/mL
NT-proANP <sub>79-98</sub> (preproANP <sub>104-123</sub> )	1.0 ng/mL
BNP <sub>32</sub> (Natreacor®)	3.5 µg/mL
CNP <sub>22</sub>	2.2 µg/mL
DNP	1.0 ng/mL
VNP	1.0 ng/mL
Adrenomedullin	1.0 ng/mL
Aldosterone	0.6 ng/mL
Angiotensin I	0.6 ng/mL
Angiotensin II	0.6 ng/mL
Angiotensin III	1.0 ng/mL
Endothelin	20 pg/mL
Renin	50 ng/mL
Urodilatin	3.5 µg/mL
Arg-Vasopressin	1.0 ng/mL

### Analytical sensitivity

The analytical sensitivity of the PBNP method is ≤10.0 pg/mL. Analytical sensitivity is defined as the concentration at two standard deviations (N=20) above a sample devoid of PBNP, such as the PBNP Level 1 Calibrator (0 pg/mL).

### Dilutional Linearity

A patient plasma sample with a high NT-proBNP concentration was diluted across the assay range. Two additional plasma samples were diluted linearly to confirm dilutional accuracy at the low end of the assay range. All samples were diluted with low level patient plasma pools. The following table summarizes the results of these studies:

	Sample A	Sample B	Sample C
Starting Concentration (pg/mL)	35979	4616	1037
<b>Regression Statistics</b>			
Slope	0.954	0.990	1.030
Intercept	511.9	-8.5	-1.6
Correlation Coefficient	0.999	1.000	0.999
N	6	6	6
<b>% Recovery</b>			
Average	99.8	98.6	101.0
Range	95.9-104.5	96.4-100.0	96.8-103.8

Manufactured under license from Roche Diagnostics.

US 5786163 and foreign equivalents

US Patent Application: 09/890442 and foreign equivalents

Symbols Key: See adjacent panel.

### Bibliography:

- Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001;22:1527-1560.
- Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998;97:1921 - 1929.
- de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230:767 - 770.
- Valli N, Gobinet A, Bordernave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999;134:437-444.
- de Bold AJ, Boerenstein HB, Varess AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981;28:89-94.
- Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987;79:738-745.
- Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999;20:1374-1375.
- Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47:287-296.
- Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999;20:1736-1744.
- Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996;78:284-287.
- McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998;351:9-13.
- Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89:879-881.
- James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108:275-281.
- Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotrophin 1 are raised in unstable angina. Heart 2000;84:421-424.
- Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002;40:437-445.
- Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004;102:37-40.
- Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003;107:2786-2792.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3<sup>rd</sup> edition, W.B Saunders Co. Philadelphia, PA 1999 pp. 42-72 (specimen collection).
- Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 982;127:308 - 315.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')<sub>2</sub> conjugate and polyclonal mouse IgG. Clin Chem 1992;38:1737-1742.
- Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.

Dimension® and Flex® are registered trademarks of Dade Behring Inc., in the U. S. Patent and Trademark Office, in Germany and many other countries.

\* An asterisk in the margin denotes a revised section.

DADE BEHRING

Dimension® clinical chemistry system  
Heterogeneous Immunoassay Module

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### Literatur:

- Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001;22:1527-1560.
- Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998;97:1921 - 1929.
- de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230:767 - 770.
- Valli N, Gobinet A, Bordernave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999;134:437-444.
- de Bold AJ, Boerstein HB, Varess AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981;28:89-94.
- Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987;79:738-745.
- Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999;20:1374-1375.
- Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47:287-296.
- Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999;20:1736-1744.
- Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996;78:284-287.
- McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998;351:9-13.
- Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89:879-881.
- James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108:275-281.
- Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotrophin 1 are raised in unstable angina. Heart 2000;84:421-424.
- Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002;40:437-445.
- Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004;102:37-40.
- Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003;107:2786-2792.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3rd edition, W.B Saunders Co. Philadelphia, PA 1999 pp. 42-72 (specimen collection).
- Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 1982;127:308 - 315.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')<sub>2</sub> conjugate and polyclonal mouse IgG. Clin Chem 1992;38:1737-1742.
- Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.

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### Bibliographie:

- Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001;22:1527-1560.
- Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998;97:1921 - 1929.
- de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230:767 - 770.
- Valli N, Gobinet A, Bordenave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999;134:437-444 .
- de Bold AJ, Boerenstein HB, Varess AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981;28:89-94.
- Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987;79:738-745.
- Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999;20:1374-1375.
- Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47:287-296.
- Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999;20:1736-1744.
- Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996;78:284-287.
- McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998;351:9-13.
- Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89:879-881.
- James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108:275-281.
- Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotropin 1 are raised in unstable angina. Heart 2000;84:421-424.
- Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002;40:437-445.
- Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004;102:37-40.
- Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003;107:2786-2792.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3rd edition, W.B Saunders Co. Philadelphia, PA 1999 pp. 42-72 (specimen collection).
- Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 982;127:308 - 315.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')<sub>2</sub> conjugate and polyclonal mouse IgG. Clin Chem 1992;38:1737-1742.
- Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.

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 Dade Behring Inc.  
Newark, DE 19714, U.S.A.

 **EC REP** Dade Behring Limited  
Regus House, Atterbury  
Milton Keynes MK10 9RG  
United Kingdom

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### Bibliografia:

- Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001;22:1527-1560.
  - Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998;97:1921 - 1929.
  - de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230:767 - 770.
  - Valli N, Gobinet A, Bordenave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999;134:437-444 .
  - de Bold AJ, Boerenstein HB, Varess AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981;28:89-94.
  - Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987;79:738-745.
  - Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999;20:1374-1375.
  - Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47:287-296.
  - Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999;20:1736-1744.
  - Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996;78:284-287.
  - McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998;351:9-13.
  - Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89:879-881.
  - James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108:275-281.
  - Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiostrophin 1 are raised in unstable angina. Heart 2000;84:421-424.
  - Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002;40:437-445.
  - Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004;102:37-40.
  - Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003;107:2786-2792.
  - Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3<sup>rd</sup> edition, W.B Saunders Co. Philadelphia, PA 1999 pp. 42-72 (specimen collection).
  - Elecsys<sup>®</sup> proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
  - Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 1982;127:308 - 315.
  - Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')<sub>2</sub> conjugate and polyclonal mouse IgG. Clin Chem 1992;38:1737-1742.
  - Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.
- Dimension<sup>®</sup> e Flex<sup>®</sup> sono un marchio registrato della Dade Behring Inc. presso l'Ufficio Brevetti e Marchi in U.S. in Germania ed in molti altri paesi.

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- Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure. Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001;22:1527-1560.
- Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998;97:1921 - 1929.
- de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230:767 - 770.
- Valli N, Gobinet A, Bordernave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999;134:437-444.
- de Bold AJ, Boerenstein HB, Vares AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981;28:89-94.
- Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987;79:738-745.
- Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999;20:1374-1375.
- Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47:287-296.
- Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999;20:1736-1744.
- Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996;78:284-287.
- McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998;351:9-13.
- Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89:879-881.
- James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108:275-281.
- Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiostrophin 1 are raised in unstable angina. Heart 2000;84:421-424.
- Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002;40:437-445.
- Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004;102:37-40.
- Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003;107:2786-2792.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 3rd edition, W.B Saunders Co. Philadelphia, PA 1999 pp. 42-72 (specimen collection).
- Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 1982;127:308 - 315.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')<sub>2</sub> conjugate and polyclonal mouse IgG. Clin Chem 1992;38:1737-1742.
- Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.

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Dimension® und Flex® sind eingetragene Warenzeichen von Dade Behring Inc., im US Patent and Trademark Office, in Deutschland sowie in vielen anderen Ländern.