Stryker France
ZAC – Avenue de Satolas Green
69330 Pusignan
t: 33 (0)4 72 45 36 00 f: 33 (0)4 72 45
36 99
prenom.nom@stryker.com

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To all user of SolisRS

Schlieren, 16th of September 2008

Prewarning letter about update of IFU of SolisRS

Dear customers

We from Degradable Solutions AG, the manufacturer of SolisRS, decided to improve the instructions for use (IFU) of the SolisRS resorbable anterior cervical cage with single use implant holder, with respect to the latest results of investigations in relation with two patient events with broken SolisRS cages involved.

Research by Degradable Solutions AG and the STEGA Foundation in Amsterdam indicates that the breakage of the cages in the two patients could have been related to very unfavorable loading situations of the implant. Both patients showed initial spinal disc height smaller than 3 mm, so that an over-distraction of the ligaments was necessary to easily introduce the implant. This resulted in a wedge shaped intervertebral space, which might have led to an unfavorable load distribution on the implant.

In the following weeks you will receive labels (stickers) containing the additional information to the instructions for use of the SolisRS. The updated version of the IFU excludes patients with spinal disc height smaller than 3 mm from the treatment with SolisRS. In the meantime please find the changes to the existing IFU in the attachment (added contra-indication, changed contra-indication, added pre-operative precaution).

This measure is a precaution taken by Degradable Solutions AG. Patients already treated with SolisRS with spinal disc height smaller than 3 mm are expected to have a successful outcome of the medical procedure, as did the two patients with broken SolisRS cages mentioned in this letter. None of the two had to undergo revision surgery.

With kind regards

Degradable Solutions AG

Dr. K. Ruffieux
Security Officer