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smiths medical

Smiths Medical ASD, Inc.
Gary, Indiana 46406, USA

Smiths Medical International, Ltd.
Hythe, Kent, CT21 6JL, UK

www.smiths-medical.com

004/009/000-002 REV. 100

2018-11

1

en	English	3
fr	Français	7
de	Deutsch	12
it	Italiano	16
es	Spanish	20
pt	Portuguese	24
nl	Nederlands	28
sv	Svenska	32
da	Dansk	36
no	Norsk	40
fi	Suomi	44
el	Ελληνικά	48
cs	Čeština	53
pl	Polski	57
hu	Magyar	62
tr	Türkçe	66

en English

Bivona® Inner Cannula for Bivona® Adult Tracheostomy Tubes

- These instructions contain important information for safe use of the product. Read the entire contents of these Instructions for Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient.
- It is the responsibility of the healthcare practitioner to assure that the instructions for use and maintenance are understood by and provided to the care giver.

1. Description

A pack of single patient use inner cannulae only for use with the below Bivona® Adult Tracheostomy tubes. The inner cannula has been determined to be MR-Safe however users should refer to the relevant Bivona® Adult Tracheostomy tube instructions for use for any MR conditions required. For further information related to the tracheostomy tube, refer to the instructions for use provided with the tracheostomy tube.

60A170, 60A180 & 60A190

Bivona® Cuffless Adult Tracheostomy Tube, sizes 7.0 – 9.0 mm

750170, 750180 & 750190

Bivona® Mid-Range Aire-Cuf® Adult Tracheostomy Tube, sizes 7.0 – 9.0 mm

670170, 670175, 670180, 670185 & 670190

Bivona® TTS™ Adult Tracheostomy Tube, sizes 7.0 – 9.0 mm

Table 1

Tracheostomy Tube Product Code	Tracheostomy Tube ID (mm)	Compatible Inner Cannula Product Codes		Inner Cannula ID (mm)	Recommended Suction Catheter Size	
		Pack of 10	Pack of 2		F	(mm)
60A170, 750170, 670170	7.0	BRC470	BRC270	5.5	10	3.33
670175	7.5	BRC475	BRC275	6.0	10	3.33
60A180, 750180, 670180	8.0	BRC480	BRC280	6.5	10	3.33
670185	8.5	BRC485	BRC285	7.0	10	3.33
60A190, 750190, 670190	9.0	BRC490	BRC290	7.5	12	4.0

2. Indications for Use

A range of sterile single patient use replacement inner cannula intended to assist in the provision of a patent airway in the tracheostomised patient when used in conjunction with Bivona® Adult cuffless, Mid-Range Aire-Cuf® and TTS™ Tracheostomy tubes.

3. Contraindications

Do not use in conjunction with lasers and electrosurgical devices where their output may contact and damage the tube.

4. Warnings

- 4.1 Ensure that the inner cannula is both the correct diameter and length for the tracheostomy tube in use. Size markings are provided on the tracheostomy tube flange, inner cannula tube, and pack labels to help identify the correct inner cannula to use. Use only an inner cannula specifically supplied for use with the Bivona® Adult tracheostomy tube range. Use of an inner cannula of the wrong diameter may cause difficulties with introduction or unnecessarily restrict gas flow. An inner cannula which is too long may protrude excessively from the outer tube, leading to tracheal damage or occlusion. Use of an inner cannula which is too short may lead to build-up of secretions within the end of the tracheostomy tube which may cause infection or occlusion.**

- 4.2 Do not fit an inner cannula when using a tracheostomy tube in conjunction with a right angled connector that incorporates a co-axial fresh gas tube which protrudes into the connector lumen (e.g. Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas, resulting in barotraumas/hypoxia.
- 4.3 Do not lubricate the inner cannula with lubricating jelly as this may occlude the inner cannula and prevent the inner cannula being retained within the tracheostomy tube.
- 4.4 The patency of the Bivona® tracheostomy tube lumen must be assured by regular suctioning and the inner cannula should be routinely checked, cleaned or replaced at regular intervals to avoid blockage or reduction in lumen of the airway.
- 4.5 If difficulty is experienced inserting the inner cannula into the tracheostomy tube, do not use excessive force as this may result in kinking of the inner cannula or trauma to the patient. The inner cannula should be discarded if damaged and a replacement inner cannula used. If the replacement inner cannula cannot be inserted, the tracheostomy tube should be replaced with a new tracheostomy tube or temporarily used without an inner cannula until other remedial actions can be taken by a physician.
- 4.6 If difficulty is experienced removing the inner cannula from the tracheostomy tube, do not use excessive force to remove it as the outer tube may be removed accidentally. Both inner cannula and tracheostomy tube should be removed together and replaced with a new tracheostomy tube and inner cannula.
- 4.7 Patient's breathing gas should be adequately humidified to minimize encrustation of the tracheostomy tube and inner cannula which can result in restricted breathing capability.
- 4.8 During and after attachment of the breathing system to the tracheostomy tube connector avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.


5. Cautions

- 5.1 Always remove the inner cannula from the tracheostomy tube before cleaning.
- 5.2 Do not use any abrasive cleaning implements to clean the inner cannula. Use sterile water, saline or cleaners listed in Section 8.3. Only use the recommended cleaning swab (Product code BICSWB). Do not use this swab to clean the tracheostomy tube. Please refer to the cleaning instructions in section 8 for further information.
- 5.3 Care should be taken to ensure that the inner cannulae do not become kinked or damaged during cleaning, and that no kinked or damaged inner cannulae are reinserted into the tracheostomy tube.
- 5.4 The security of all breathing system connectors should be checked when the circuit is established and frequently thereafter.
- 5.5 Replacement inner cannula for Bivona® Adult Tracheostomy Tubes are designed and intended for single patient use. The manufacturer recommends that the inner cannula usage does not exceed twenty-nine (29) days. They should be cleaned in accordance with the instructions provided and should not be re-sterilised or cleaned by any other method (e.g. autoclave or dishwasher) and should not be placed in another patient.
- 5.6 Federal (USA) law restricts this device to sale by or on the order of a physician.
- 5.7 Dispose of this product in a safe manner according to local guidelines for disposal of contaminated waste.
- 5.8 Follow universal precautions as specified by the Centers for Disease Control and Prevention, USA, or relevant advisory bodies in the country where this device is being used.
- 5.9 When using the tracheostomy tube with inner cannula there is a reduction in the lumen size which could increase the pressure drop. If concerned, consult your healthcare provider regarding whether a larger sized tube is appropriate.

- 5.10 If suctioning the inner cannula ensure the suction catheter easily fits within the internal diameter of the inner cannula. Table 1 provides recommended suction catheter sizes to use but these are given as guidance only and the advice of clinical professionals should be followed.

6. Instructions for Use

- 6.1 Check the inner cannula prior to insertion/reinsertion and discard if kinked or damaged.
- 6.2 Insert the correct size inner cannula into the tracheostomy tube. During insertion of the inner cannula align the hub of the inner cannula to match the internal profile of the tracheostomy tube connector. Be sure that the inner cannula is securely inserted. A slight snap should be felt indicating proper positioning, if no snap is felt discard the inner cannula and replace with a new inner cannula. The inner cannula hub is designed to work with the connector design labeled (a) in the image below. Previous connector designs such as that labeled (b) in the image below will not accept the inner cannula fully.

7. Replacing the Inner Cannula

- 7.1 Disconnect the tracheostomy tube from the breathing system if attached.
- 7.2 Remove the existing inner cannula by means of the grey ring pull on the inner cannula.
- 7.3 Insert a new or cleaned inner cannula. A slight snap should be felt indicating proper positioning, if no snap is felt discard the inner cannula and replace with a new inner cannula.
- 7.4 Reconnect the breathing system, if necessary, and check for security.
- 7.5 Discard the old inner cannula if damaged.

8. Cleaning the Inner Cannula

Daily cleaning of the inner cannula is recommended or whenever the inner cannula becomes contaminated by mucus or secretions. This period of time will vary according to the individual patient's needs. Patients who have thick, tenacious, or large volumes of secretions may require more frequent inner cannula replacements.

Bivona® Inner Cannulae are thin walled and constructed of a soft and flexible material (ePTFE) and should be handled and cleaned with care. If difficulty is experienced when inserting the inner cannula into the tracheostomy tube, do not use excessive force as this may cause compression of the inner cannula. If this occurs, discard the damaged inner cannula and use a new one. Use of an inner cannula which is too short may result in secretions being left in the tube even when the inner cannula is removed. This may lead to a build-up of secretions within the end of the tracheostomy tube which may cause infection or occlusion. Do not re-sterilize or use brushes, autoclaves, dishwashers, vigorous or aggressive cleaning methods.

Only use the provided cleaning swab. Do not use any other type of swab or brush.

The following methods must be followed:

- 8.1 To clean the Bivona inner cannula firstly remove it from the tracheostomy tube by holding the grey ring pull and withdrawing it gently away from the tube.
- 8.2 Rinse with water.

- 8.3 Soak the inner cannula in a non-coloured, non-scented enzymatic cleaner (e.g. Ruhof 345PANS Endozime AW No Scent) to remove all biological material. Refer to the cleaning solution manufacturer's recommendations for soaking times, use and handling of the cleaner.
- 8.4 If enzymatic cleaner is not available, then soaking in warm water with mild soap solution may soften secretions to a level at which they can be removed by gentle swabbing.
- 8.5 Rinse the inner cannula well with clean water.
- 8.6 Following rinsing, insert the dedicated cleaning swab into the inner cannula at the ring pull end and use a twisting motion as you advance the swab to remove any contamination.
- 8.7 Continue to push the cleaning swab through the far end of the inner cannula in a forward twisting motion. Do not pull it back through the inner cannula.
- 8.8 If, following the instructions above, the inner cannula is still contaminated, then the secretions cannot be removed without unacceptable risk of damaging the inner cannula and it should be discarded and a new inner cannula used.
- 8.9 Dispose of cleaning swab in a safe manner according to local hospital procedures.
- 8.10 Rinse inner cannula with water.
- 8.11 Allow the inner cannula to air dry. Once clean and dry it should be stored free from particulate matter.
- 8.12 If any dust, contamination, damage or kinking to the inner cannula is observed before replacing the inner cannula, discard the inner cannula and replace it with a new one.

The Smiths Medical and Portex design marks; Aire Cuf; TTS; and Bivona are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Pat pending.

Consult instructions for use • **Caution** • Do not re-sterilize • Does not contain Natural Rubber Latex • Do not use if package is damaged • Sterilized using ethylene oxide • MR Safe • Catalogue Number • Batch Code • Date of Manufacture • Use by • Manufacturer • Quantity

© 2011, 2018 Smiths Medical family of companies. All rights reserved.

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Smiths Medical og Portex designmærkerne, Aire Cuf, TTS samt Bivona er varemærker, der tilhører Smiths Medical virksomhedsgruppen. Symbolet ® angiver, at varemærket er registreret hos det amerikanske patent- og varemærkekontor (U.S. Patent and Trademark Office) og i visse andre lande. Patentanmeldt.
Znaki graficzne Smiths Medical i Portex; Aire Cuf; TTS; i Bivona są znakami towarowymi przedsiębiorstw Smiths Medical. Symbol ® wskazuje, że znak towarowy jest zarejestrowany w amerykańskim urzędzie ds. patentów i znaków towarowych (U.S. Patent and Trademark Office) oraz w niektórych innych krajach. Patenty w toku.