# 3PLUS1<sup>®</sup> Multilevel Plasma Calibrator Set MassTox<sup>®</sup> Neuroleptics 1/*EXTENDED*

Order no. 92028/XT

## Instructions

These Chromsystems plasma calibrators are designed for calibrating your analysis system. The lyophilised calibrators are based on human plasma and should be used routinely according to the assay procedure being used. They are available as Multilevel Calibrator Set only, with 3 different calibration levels and a blank. The calibrators are to be handled in the same manner as a patient specimen.

## Reconstitution

Carefully remove the rubber stopper from the vial. Add exactly 1.0 ml of distilled water into the vial. Replace the stopper and allow to sit for 10 to 15 minutes at room temperature. Swirl the vial to dissolve the contents until homogeneity.

## Storage conditions

Stored below -18 °C, unopened calibration standards are stable until the expiry date which is stated on the vial labels.

Light protected and tightly capped, the reconstituted calibrators are stable for 14 days at +2 to +8 °C except Olanzapine\*\* and Desmethylolanzapine\*\* which should be aliquoted immediately after reconstitution and stored below -18 °C. If the calibrators cannot be analysed within this period, they should be aliquoted and stored frozen below -18 °C (up to 3 months).

## Caution

This product is manufactured from pooled human plasma which has been tested by the manufacturer and found negative for infections by the human immunodeficiency virus (HIV), the hepatitis B virus (HBV), the hepatitis C virus (HCV) and the bacterium *Treponema pallidum*. Nevertheless, a potential risk of infection cannot be entirely excluded. Consider all products containing human source material as potentially infectious and exercise the same care in the handling of this product as in the handling of potentially infectious patient samples.

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**\*\*Note:**

These analytes are not stable at +2 to +8 °C and at room temperature. The calibrators should be aliquoted immediately after reconstitution and stored below -18 °C.

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**Traceability/Traçabilité**

The values of this calibrator of Dehydroaripiprazole, Clozapine, N-Desmethylclozapine, Haloperidol, Olanzapine, Quetiapine, Norquetiapine, Risperidone and 9-OH-Risperidone are traceable to certified reference materials and primary standards.

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