Rév 1 : Mars 2020

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Réf ACST : 455 - CAPA N1

Date : 16 Mars 2020

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INSTRUCTIONS FOR USE

DESCRIPTION

Vantris VUR is a non-absorbable, permanent-action, tissue bulking substance, for the endoscopic treatment of vesicoureteral reflux (VUR). It is easily and manually injected.

Vantris VUR consists of particles of polyacrylate polyalcohol immersed in a glycerol and saline solution carrier. It has a very high molecular mass (~10 million Daltons) and it comes in the form of sterile, pyrogen-free particles that are highly deformable by compression, and that may be injected with 23-gauge needles.

Once implanted, no local, regional or distance migration has been observed.

The carrier is a 40% glycerol solution. Once implanted, it is eliminated by the reticuloendothelial system without metabolizing and excreted through the kidneys, while the particles remain to ensure permanent bulking.

INDICATIONS

The physicochemical properties of Vantris VUR make it suitable for the endoscopic treatment of vesicoureteral reflux (VUR).

MECHANISM OF ACTION

When Vantris VUR is implanted in the ureterovesical junction, the material enlarges the volume of the area and corrects the anatomy of the meatus and the distal ureter, preventing urine to return to the ureter after having been stored in the bladder.

HOW SUPPLIED

Vantris VUR (code BAR- 1J) is supplied as follows:

- 1 syringe (Code: BARI-1J) with stopper, containing 1 ml of material. Vantris VUR is provided sterile and ready-to-use. The syringe is sterilized with moist heat.

Vantris VUR can be injected using either a flexible injection needle or a metal semi-rigid injection needle.

It is recommended to use the injection needles developed by Promedon, with the following features:

- 1 23 G flexible injection needle (Code 50F or 37F)

Or alternatively:

- 1 22 G metal semi-rigid injection needle (Code RIN or RINS)

PACKAGING

Vantris VUR is supplied in a box containing the syringe in a triple pouch; the external one is made of aluminum.

After use, discard the product and packaging according to hospital, administrative or local government procedure.

HANDLING AND STORAGE PRECAUTIONS

Vantris VUR have been designed to be used only ONCE. Therefore, DO NOT REUSE or RESTERILIZE, as this can potentially result in compromised device performance and increased risk of inappropriate resterilization and cross contamination.

IF ANY POUCH CONTAINING Vantris VUR COMPONENTS WERE DAMAGED, DO NOT IMPLANT.

IMPORTANT: Only the internal pouch can enter the sterile area.

STORAGE

Vantris VUR should be stored under the following conditions:

- TEMPERATURE: Room temperature.

DO NOT USE AFTER THE EXPIRATION DATE ON THE PACKAGING.

OPERATING ROOM CONDITIONS

Operating room conditions must meet hospital, administrative or local government procedure.

PATIENT INFORMATION

Surgeons are responsible for informing patients or their representatives on the possible complications associated to the Vantris VUR implant surgery, in advance.

CONTRAINDICATIONS

- Ureterocele
- Voiding dysfunction
- Kidney dysfunction

- Paraureteral (Hutch) diverticulum
- Urinary tract infections
- Ureterovesical junction obstruction

PRECAUTIONS

The technique and the instruments necessary for injecting Vantris VUR should follow the same precautions as in other similar urological interventions, endoscopy practices and endourological procedures.

Vantris VUR should be used by qualified surgeons with experience in endoscopy practices and trained in subureteral injection procedures.

The administration of prophylactic therapy with antibiotics should be considered, according to the procedure approved by the hospital.

WARNINGS

Tissues with exacerbated scarring processes and/or with fibrosis can compromise the ideal bulking.

Vantris VUR should not be injected intravascularly since it may cause vascular occlusion.

The maximum amount of Vantris VUR Treatment product to be implanted in a single surgery of Renal Refluxing Unit (RRU) must not exceed 1 KIT.

POSSIBLE COMPLICATIONS

The possible complications associated with the use of Vantris VUR should be discussed with the patient's parents prior to the surgery.

Complications associated with the implantation of an injectable agent into the vesicoureteral junction include those common to this types of procedures:

- Temporary post-injection dilatation of the upper urinary tract.
- Obstruction of the vesicoureteral junction.
- Dysuria
- Hematuria
- Infection
- Technique failure

Promedon requires that all surgeons inform the Company of any complication observed with the use of Vantris VUR.

MAGNETIC RESONANCE (MR) ENVIRONMENT

The implant does not affect and is not affected by magnetic resonance (MR) environments.

IMPLANTATION PROCEDURE

In order to perform the procedure, the following components are necessary:

- 1 syringe (Code: BARI-1J): The quantity to be used will depend on each particular case.
- 1 flexible injection needle or alternatively 1 metal semi-rigid injection needle.

It is recommended to use the injection needles developed by Promedon, with the following features:

- 1 23 G flexible injection needle (Code 50F or 37F)

Or alternatively:

- 1 22 G metal semi-rigid injection needle (Code RIN or RINS)
- 1 cystoscope with a straight working channel of 4 Fr or more, according to the needle chosen.

1. Vantris VUR is injected under general anesthesia in children.
2. The free flow through the injection needle is verified using saline solution.
3. The syringe is connected to the injection needle passing material through it until the substance appears at the needle tip. Insert the injection needle into the cystoscope.
4. Prior to the injection, the meatus should be observed with different volumes of bladder filling in order to choose the ideal situation.
5. The submucosa of the bladder is punctured at the six o'clock position, 3 mm under the ureteral meatus, 4 to 5 mm depth.
6. The material is then slowly injected until the ureteral wall is adequately modified (bulkiness). More punctures may be performed until the desired effect is reached.
7. Once the injection is completed, the needle is kept in its position for 30 seconds.
Then, the needle is removed.
8. Once the procedure has been completed, the bladder is emptied, the optic removed, and an optional voiding cystoureterogram (VCUG) may be performed to verify reflux repair.



NOTES

Consider that 0.4 ml of the product remains in the injection needle, which can be used by pushing it down with the same syringe with sterile physiological solution or sterile distilled water.

A ureteral catheter may be necessary where the meatus presents difficulties for injection.

Re-implantation needs: If after performing a voiding cystoureterogram (VCUG) three months after the treatment the persistence of reflux is proven, Vantris VUR can be re-injected until it is eliminated.

SYMBOLS USED IN LABELS

	PRODUCT CODE
	LOT NUMBER
	CAUTION
	EXPIRY DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE. STERILIZATION METHOD: MOIST HEAT
	DO NOT RESTERILIZE
	MANUFACTURER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	SERIAL NUMBER
	READ THE INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGE IS DAMAGED
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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