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Dimension® clinical chemistry system

Flex® reagent cartridge

IBCT

See shaded sections: Updated information from 2008-03 version.

Issue Date 2009-10-07

Total Iron Binding Capacity (IBCT)


Intended Use: The IBCT method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure Total Iron Binding Capacity in human serum.

Summary: Total iron binding capacity (IBCT) is a measure of the serum transferrin iron binding capacity. Measurements of serum iron and total iron binding capacity are widely used in the diagnosis and treatment of iron deficiency anemia and chronic inflammatory disorders.¹

Principles of Procedure: This fully automated method adds saturating iron to the sample to saturate the transferrin iron-binding sites. The excess unbound iron is photometrically determined (instead of being physically removed by adsorption), in a manner similar to those described by Yamanishi et. al.² Subsequent addition of acid causes the release of bound iron from transferrin, which is then analyzed using the chromogen Ferene®. Earlier work by Higgins,³ Artiss, et. al.⁴ and Hennessy et. al.⁵ demonstrated the high sensitivity of Ferene® and its utility in iron assays. A surfactant is used to prevent protein precipitation. Reaction times are optimized to minimize interferences, and thiourea is present to complex copper.

The serum sample is automatically mixed with a ferric iron solution, which saturates all available iron-binding sites of transferrin. Under non-acidic conditions (pH 8.6), only unbound, excess saturating iron is available to be reduced to ferrous iron by ascorbic acid and to form a blue complex with Ferene®. Subsequent addition of acid (final pH of 4.5) releases the iron bound to transferrin; this additional iron is reduced to ferrous iron by ascorbic acid and forms an increased amount of blue complex with Ferene®. The increase in absorbance upon shifting from pH 8.6 to pH 4.5, measured using a bichromatic (600, 700 nm) endpoint technique, is proportional to the concentration of transferrin-bound iron, and thus to the iron binding capacity (total) of the serum sample.

Ferene® is a registered trademark of Diagnostic Chemicals, LTD., Charlottetown, P.E.I., Canada C1A4H5, for the compound 3-(2-pyridyl) 5,6-bis-2-(5-furyl sulfonic acid)-1, 2,4-triazine disodium salt.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b
1, 2, 3	Tablet	Ascorbic acid	19 mM
4	Liquid	Ferene®	0.56 mM
5, 6	Liquid	Ferric chloride	0.02 mM
		Citric acid	0.2 mM
7	Liquid	Acetate Buffer	500 mM
		Thiourea	33 mM
8	Liquid	Tris buffer	200 mM

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value in final reaction mixture.

Risk and Safety:

Harmful. Contains thiourea.
Contains mixture of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

- R40: Limited evidence of a carcinogenic effect
R36/38: Irritating to eyes and skin.
R43: May cause sensitisation by skin contact.
S24: Avoid contact with skin.
S37: Wear suitable gloves.

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the Dimension® system.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 6 days for wells 1 – 3
30 days for wells 4 – 8

Specimen Collection and Handling: Normal procedures for collecting serum may be used for samples to be analyzed by this method.⁶

Follow the instructions provided with your specimen collection device for use and processing.⁷

For serum, complete clot formation should take place before centrifugation.⁸

Serum should be removed from cells within two hours.⁹ Separated specimens are stable for 4 days at room temperature and up to 7 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C for up to 2 months.¹⁰

Procedure

Materials Provided

IBCT Flex® reagent cartridge, Cat. No. DF84

Materials Required But Not Provided

IBCT Calibrator, Cat. No. DC84
Quality Control Materials

Test Steps

Sampling, reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

Test Conditions

Sample Size	25 µL
Reagent 1 Volume	36 µL
Reagent 2 Volume	25 µL
Reagent 3 Volume	25 µL
Reagent 4 Volume	75 µL
Temperature	37 °C ± 0.1 °C
Wavelength	600 and 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	0 – 1000 µg/dL [0 – 179.1 µmol/L] ^c
Calibration Material	IBCT Calibrator, Cat. No. DC84
Calibration Scheme	3 levels, n = 3
Units	µg/dL [µmol/L] (µg/dL x 0.179) = [µmol/L]
Typical Calibration Levels	0, 550, 1100 µg/dL [0, 98, 197 µmol/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	<ul style="list-style-type: none"> For each new lot of Flex® reagent cartridges After major maintenance or service, if indicated by quality control results As indicated in laboratory quality control procedures When required by government regulations
Assigned Coefficients	C ₀ -50.00 C ₁ 7.0

c. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known total iron binding capacity concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of Total Iron Binding Capacity (IBCT) in µg/dL [µmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0 – 1000 µg/dL [0 – 179.1 µmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 1000 µg/dL [179.1 µmol/L] should be repeated on dilution.

Manual Dilution: Make appropriate dilution with Reagent grade water to obtain results within reportable range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide. Recommended Auto Dilute volume is 12 µL.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
550 µg/dL [98.5 µmol/L]	>9.5 µg/dL [1.7 µmol/L]
1000 µg/dL [179.1 µmol/L]	>22.1 µg/dL [3.96 µmol/L]

Interfering Substances

Measurements of IBCT may be inaccurate if performed within 14 days of IV iron dextran administration.¹¹ Ferrus sulfate at 250 µg/dL increased IBCT results by 17%.

Hemoglobin (hemolysate) of 300 mg/dL [1.86 mmol/L] increased IBCT results at 388 µg/dL [69 µmol/L] by 14%.

Expected Values: 250 – 450 µg/dL [44.8 – 80.6 µmol/L]¹²

Each laboratory should establish its own reference interval for total iron binding capacity as performed on the Dimension® system.

Specific Performance Characteristics^d

Precision^{e,f}

Material	Mean µg/dL [µmol/L]	Standard Deviation (% CV)	
		Within-run	Total
Serum Pool 1	192 [34.4]	4.1 [0.7] (2.1)	6.2 [1.1] (3.2)
Serum Pool 2	329 [58.9]	5.4 [1.0] (1.6)	9.5 [1.7] (2.9)
Serum Pool 3	513 [91.9]	5.9 [1.1] (1.2)	12.2 [2.2] (2.4)
Liquichek™ L 1	328 [58.7]	5.5 [1.0] (1.7)	8.9 [1.6] (2.7)
Liquichek™ L 2	206 [36.9]	4.1 [0.7] (2.0)	6.9 [1.2] (3.3)

d. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (Refer to your Dimension® Operator's Guide).

e. Reproducibility testing was done in accordance with the NCCLS/CLSI Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999).

f. Specimens at each level were analyzed in duplicate, twice a day for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

Liquichek™ is a trademark of Bio-Rad Diagnostics Group, Irvine, CA 92714.

Method Comparison

Regression Statistics^g

Comparative Method	Slope	Intercept µg/dL [µmol/L]	Correlation Coefficient	n ^h
Dimension® TIBC	1.14	-22.4 [-4.0]	0.971	137

g. Model equation for regression statistics is: [result of IBCT] = slope x (comparative method result) + intercept.

h. The range of TIBC values in the correlation study was 59 to 469 µg/dL [10.6 to 84.0 µmol/L].

Specificity

HIL Interference

The IBCT method was evaluated for interference from hemolysis, icterus, and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	IBCT Concentration µg/dL [µmol/L]	Bias ⁱ %
Hemoglobin (hemolysate)	200 mg/dL [0.31 mmol/L] (monomer)	388 [69]	<10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	386 [69]	<10
Lipemia (Intralipid®)	600 mg/dL [6.78 mmol/L]	383 [69]	<10

i. Analyte results should not be corrected based on this bias.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Non-Interfering Substances

The following substances do not interfere with the IBCT method when present in serum at the concentrations indicated. Systemic inaccuracies (bias) due to these substances are less than 5% at a total iron binding capacity level of 260 µg/dL [47 µmol/L]:

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic Acid	3 mg/dL	170.3 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	20 ng/mL	25.6 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	800 mg/dL	174 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Ferritin	12000 ng/mL	12000 µg/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	8 U/mL	8000 U/L
Ibuprofen	40 mg/dL	1939 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium	3.5 mg/dL	5.04 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Protein, human albumin	6 g/dL	60 g/L
Protein, human IgG	6 g/dL	60 g/L
Salicylic acid	100 mg/dL	7.24 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 6 µg/dL [1.1 µmol/L]

The analytical sensitivity represents the lowest concentration of total iron binding capacity that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics Inc.
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics Inc.
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía

1. Tietz NW. Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Burtis CA, Ashwood ER, Bruns DE, Elsevier Saunders, 2006; pp 1186-1193.
2. Yamanishi H. et al. Fully automated measurement of total iron-binding capacity in serum. Clin Chem 43:12, 2413-2417;1997.
3. Higgins T. Novel chromogen for serum iron determinations. Clin Chem 1981;27:1619-1620.
4. Artiss JD, Vinogradov S, Zak B. Spectrophotometric study of several sensitive reagents for serum iron, Clin Biochem 1981;14:311-315.
5. Hennessy DG, Reid GR, Smith FE, Thompson SL. Ferene—a new spectrophotometric reagent for iron. Can J Chem 1984;62:721-724.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
7. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
8. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
9. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1999, Third Edition pp 54-56.
10. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007;p 570.
11. Clinical Practice Guidelines for the Treatment of Anemia of Chronic Renal Failure; Guideline 5. National Kidney Foundation, Inc., 1997.
12. Tietz NW. Clinical Guide to Laboratory Tests W.B.Saunders Co., Philadelphia, PA 1990.

Symbols Key Symbolschlüssel Explication des Symboles Interpretazione simboli Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
EXP CCYY-MM-DD	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbricante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente para "n" saggi / Contenido suficiente para <n> ensayos
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
CONTENTS	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel

2006-06_EFIGS