# Notification urgente

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**Date** 26 mars 2010

**Produit**

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## Annexe A : Paramètres ARCHITECT cSystems (Conv / SI)

Paramètre	Valeur précédente	Nouvelle valeur	Exemple																																					
Echantillon (STANDARD)	4.0	2.6	<div style="border: 1px solid black; padding: 5px;"> <p style="margin: 0;"><b>Configure assay parameters – General</b></p> <p style="margin: 0;"><input type="radio"/> Reaction definition    <input checked="" type="radio"/> Reagent / Sample    <input type="radio"/> Validity checks</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">Reagent: <b>BILIT</b></td> <td style="width: 25%;">Reagent volume: <b>160</b></td> <td style="width: 25%;">R1 R2</td> </tr> <tr> <td>Diluent: <b>Saline</b></td> <td>Water volume: _____</td> <td>_____</td> </tr> <tr> <td>Diluent dispense mode: <b>Type 0</b></td> <td>Dispense mode: <b>Type 0</b></td> <td><b>Type 0</b>    <b>Type 0</b></td> </tr> </table> <table style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr> <th style="width: 15%;">Dilution name</th> <th style="width: 15%;">Sample</th> <th style="width: 15%;">Diluted sample</th> <th style="width: 15%;">Diluent</th> <th style="width: 15%;">Water</th> <th style="width: 15%;">Dilution factor</th> <th style="width: 15%;">Default dilution</th> </tr> </thead> <tbody> <tr> <td><b>STANDARD :</b></td> <td><b>2.6</b></td> <td>_____</td> <td>_____</td> <td>_____</td> <td><b>1:1.00</b></td> <td><input checked="" type="radio"/></td> </tr> <tr> <td>1:5 :</td> <td>20.0</td> <td><b>2.6</b></td> <td><b>80</b></td> <td>_____</td> <td><b>1:5.00</b></td> <td><input type="radio"/></td> </tr> <tr> <td>1:10 :</td> <td>10.0</td> <td><b>2.6</b></td> <td><b>90</b></td> <td>_____</td> <td><b>1:10.00</b></td> <td><input type="radio"/></td> </tr> </tbody> </table> <p style="margin: 5px 0 0 0;"><input type="radio"/> Reaction definition    <input type="radio"/> Reagent / Sample    <input checked="" type="radio"/> Validity checks</p> <p style="margin: 0 0 0 20px;">Reaction check:    <b>End Subtraction</b></p> <p style="margin: 5px 0 0 100px;">Read time:    <b>A</b>    <b>B</b></p> <p style="margin: 0 0 0 100px;">Calculation limits:    <b>-0.1000 – 0.0045</b></p> </div>	Reagent: <b>BILIT</b>	Reagent volume: <b>160</b>	R1 R2	Diluent: <b>Saline</b>	Water volume: _____	_____	Diluent dispense mode: <b>Type 0</b>	Dispense mode: <b>Type 0</b>	<b>Type 0</b> <b>Type 0</b>	Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution	<b>STANDARD :</b>	<b>2.6</b>	_____	_____	_____	<b>1:1.00</b>	<input checked="" type="radio"/>	1:5 :	20.0	<b>2.6</b>	<b>80</b>	_____	<b>1:5.00</b>	<input type="radio"/>	1:10 :	10.0	<b>2.6</b>	<b>90</b>	_____	<b>1:10.00</b>	<input type="radio"/>
Reagent: <b>BILIT</b>	Reagent volume: <b>160</b>	R1 R2																																						
Diluent: <b>Saline</b>	Water volume: _____	_____																																						
Diluent dispense mode: <b>Type 0</b>	Dispense mode: <b>Type 0</b>	<b>Type 0</b> <b>Type 0</b>																																						
Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution																																		
<b>STANDARD :</b>	<b>2.6</b>	_____	_____	_____	<b>1:1.00</b>	<input checked="" type="radio"/>																																		
1:5 :	20.0	<b>2.6</b>	<b>80</b>	_____	<b>1:5.00</b>	<input type="radio"/>																																		
1:10 :	10.0	<b>2.6</b>	<b>90</b>	_____	<b>1:10.00</b>	<input type="radio"/>																																		
Echantillon dilué (1:5)	4.0	2.6																																						
Echantillon dilué (1:10)	4.0	2.6																																						
Contrôles de validité / Limites de calculs	0.0070	0.0045																																						
Echantillon (Blanc)	4.0	2.6	<div style="border: 1px solid black; padding: 5px;"> <p style="margin: 0;"><b>Configure assay parameters – Calibration</b></p> <p style="margin: 0;"><input type="radio"/> Calibrators    <input checked="" type="radio"/> Volumes    <input type="radio"/> Intervals    <input type="radio"/> Validity checks</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">Calibrator: <b>Bil</b></td> <td style="width: 15%;">Calibrator level</td> <td style="width: 15%;">Sample</td> <td style="width: 15%;">Diluted sample</td> <td style="width: 15%;">Diluent</td> <td style="width: 15%;">Water</td> </tr> <tr> <td>Blank: <b>Water</b></td> <td></td> <td><b>2.6</b></td> <td>_____</td> <td>_____</td> <td>_____</td> </tr> <tr> <td>Cal 1: <b>Bil1</b></td> <td></td> <td><b>2.6</b></td> <td>_____</td> <td>_____</td> <td>_____</td> </tr> <tr> <td>Cal 2: <b>Bil2</b></td> <td></td> <td><b>2.6</b></td> <td>_____</td> <td>_____</td> <td>_____</td> </tr> </table> </div>	Calibrator: <b>Bil</b>	Calibrator level	Sample	Diluted sample	Diluent	Water	Blank: <b>Water</b>		<b>2.6</b>	_____	_____	_____	Cal 1: <b>Bil1</b>		<b>2.6</b>	_____	_____	_____	Cal 2: <b>Bil2</b>		<b>2.6</b>	_____	_____	_____													
Calibrator: <b>Bil</b>	Calibrator level	Sample		Diluted sample	Diluent	Water																																		
Blank: <b>Water</b>		<b>2.6</b>		_____	_____	_____																																		
Cal 1: <b>Bil1</b>		<b>2.6</b>	_____	_____	_____																																			
Cal 2: <b>Bil2</b>		<b>2.6</b>	_____	_____	_____																																			
Echantillon (Cal 1)	4.0	2.6																																						
Echantillon (Cal 2)	4.0	2.6																																						