Urgent NOTIFICATION DE SECURITE

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CAUTION notice for the Manual of the CDI 100 System

CDI™ 100 Hematocrit/Oxygen Saturation Monitoring system

Terumo CDI™ Hematocrit/Oxygen saturation monitoring systems are intended to monitor blood values including Oxygen Saturation, Hematocrit (Hct), and Hemoglobin (Hgb). When used in accordance with their instructions for use, the systems have been demonstrated to provide reliable reports of these values with an accuracy characterized in Appendix A of the Operator's Manual.

>Caution

Failure to follow the instructions can cause the monitoring system to display inaccurate values.

The accuracy of the results is dependent upon the following:

- Reading and understanding the instructions for use
- Proper set-up of the system
- Use of all available system features
- Periodic comparison to a laboratory reference sample

When a displayed value is significantly different from expectations based on the clinical situation, verify its accuracy by independent means before initiating treatment.

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*CAUTION notice for the Manual of the **CDI 101** System revisions **A to C**.*

CDI™ 101 Hematocrit/Oxygen Saturation Monitoring system

Terumo CDI™ Hematocrit/Oxygen saturation monitoring systems are intended to monitor blood values including Oxygen Saturation, Hematocrit (Hct), and Hemoglobin (Hgb). When used in accordance with their instructions for use, the systems have been demonstrated to provide reliable reports of these values with an accuracy characterized in Appendix B of the Operator's Manual.

>Caution

Failure to follow the instructions can cause the monitoring system to display inaccurate values.

The accuracy of the results is dependent upon the following:

- Reading and understanding the instructions for use
- Proper set-up of the system
- Use of all available system features
- Periodic comparison to a laboratory reference sample

When a displayed value is significantly different from expectations based on the clinical situation, verify its accuracy by independent means before initiating treatment.

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