---

MERIT MEDICAL SYSTEMS, INC  
1600 WEST MERIT PARKWAY  
SOUTH JORDAN, UTAH 84095  
PHONE 801-253-1600  
FAX 801-253-1688  
www.merit.com

**URGENT FIELD SAFETY NOTICE (FSN)**

September 8, 2016

ATTENTION: Risk Manager/Cardiac Cath Lab Manager

**Product:** Drip Chambers

**Affected Products:** See attached table.

Dear Sir or Madame,

Merit Medical Systems, Inc. is voluntarily conducting a field correction due to a supplier molding defect that has been identified in the drip chamber component of certain Fluid Administration Sets.

**This letter updates the letter you previously received by identifying the risk associated with continued use of a faulty drip chamber after the defective device has been identified (see bold text below). There have been no other changes to this letter.**

The molding defect occurs in the fluid administration spike between the vent channel and the drip channel. When the fluid administration set is prepped and contrast is drawn out of the contrast bottle, air may be drawn into the contrast fluid channel from the vent channel resulting in the contrast having excessive air bubbles, which can cause the burette to empty quickly and the blue ball to drop to the bottom of the chamber interfering with contrast administration. Merit has identified 822 lots and 130 catalog numbers (107,331 units), as shown in the attached table. **Although the defect is easily and readily identifiable during standard clinical preparation, continued use after identifying a faulty drip chamber may result in an air embolus.** Merit has received no reports of patient harm or injury as a result of this issue. Merit is providing labels and label application instructions that direct the user to evaluate the drip chamber at the point of use and, if necessary, replace the impacted fluid line including this component with a replacement product (K08-02385) that is being provided separately at no charge.

Our records indicate that you have received affected Drip Chamber lots. This letter is being provided to the Risk Manager and Cardiac Cath Lab Manager. Only the Cardiac Cath Lab Manager will receive the Customer Response Form, Product Labels, and Label Application Instructions.

Actions required of the Cardiac Cath Lab Manager:

1. Please determine if any of the devices are within your facility using the enclosed table of affected drip chamber lots and ensure that all personnel to whom the devices were distributed are made aware of this field action.
2. If any of the affected devices are within your facility, please place one of the enclosed labels next to the current pouch label of each device according to the enclosed instructions. The label directs the user on evaluating the drip chamber at point of use.
3. Complete the enclosed Customer Response Form and indicate whether or not you have any of the affected devices in your facility and the actions completed to label the affected devices.
4. Email the completed Customer Response Form to [RESPONSE@merit.com](mailto:RESPONSE@merit.com) or fax to +31433588221.
5. If you have any questions concerning this communication or for assistance labeling your affected products, please don't hesitate to contact your Merit Sales Representative or Layla Jitou at [ventes@merit.com](mailto:ventes@merit.com) or at 0800 91 60 30.

Merit Medical is committed to providing high quality products to you and apologizes for any inconvenience this field action may cause.

Kind Regards,

Mark Mullaney  
Director, Regulatory Affairs  
Merit Medical Ireland Ltd.

Enclosures: Table of Affected Drip Chamber Lots

Only included with Cardiac Cath Lab Manager Letter  
Customer Response Form  
Product Labels  
Label Application Instructions

MERIT MEDICAL SYSTEMS, INC  
1600 WEST MERIT PARKWAY  
SOUTH JORDAN, UTAH 84095  
TÉL. : 801-253-1600  
FAX : 801-253-1688  
www.merit.com

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