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Contacts:  
Pascal Bernard pbernard@ioltech.com  
Tel. +33 (0)5 46 44 08

RE : Definitive stop of marketing, product recall and follow-up of implanted patients  
Presbyopic intraocular lenses NEWLIFE® / VIVARTE® PRESBYOPIC

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**ACTION REQUIRED**

Dear Doctor,

This message is a very important update from our March 2006 communication and will require action on your part.

**Concerning presbyopic intraocular lenses, NEWLIFE® / VIVARTE® PRESBYOPIC**

As background, in March we informed you of concerns brought to our attention regarding endothelial cell loss with the Newlife® or Vivarte® Presbyopic intraocular lens, both implants being presbyopic lenses implanted in the anterior chamber. At that time we reported endothelial cell loss that started occurring after two to three years post-implantation. In these patients the loss was great enough that it led to explantation of the lens. Consequently to these incidents, Ioltech discontinued the sales and marketing of the Newlife® intraocular lens (marketing of Vivarte® Presbyopic intraocular lens was stopped in December 2005).

Since March, in agreement with our local French authority, the French Health Products Safety Agency (AFSSAPS), a retrospective investigation involving 150 patients (270 eyes) having received at least one implant Newlife® or Vivarte® presbyopic with a mean follow-up of 18.76 months [0.26-54.47] has been performed. This study shows that the risk of endothelial cell loss increases after two to three years of implantation. The level of cell loss may necessitate explantation of the lens. However, no specific cause for this endothelium cell loss has been clearly identified.

The study results support our decision (see attached synopsis). As a result, sales and marketing of Newlife® will not be resumed again.

**Concerning myopic intraocular lenses, GBR® / VIVARTE® MYOPIC**

Two other implants used in myopic correction, GBR® and Vivarte® myopic, have a similar design and are also implanted in the anterior chamber in phakic eyes. Several cases of endothelial cell loss that may necessitate explantation, although less frequent, have been reported so far.

In agreement with our local French authority AFSSAPS, Ioltech voluntarily and definitively stops selling the GBR® (marketing of Vivarte® myopic intraocular lens was stopped in December 2005).

One of our staff will be contacting your practice to determine if you have any remaining product still in stock and to facilitate the return to Ioltech as soon as possible.

**Follow-up of implanted patients**

In agreement with our local French authority AFSSAPS, for the patients having been implanted so far, we recommend that the periodic monitoring of corneal endothelium performed by cell count using specular microscopy after implantation is increased from yearly to every 6 months. In the case that you observe greater than a 30% endothelial cell loss (when compared to the preoperative cell count) and/or count below 1500 cells/mm<sup>2</sup>, the opportunity of an explantation should be evaluated in concert with the patient ; explantation may be uni- or bilateral.

Any loss of endothelial cells at a level such as described here above and/or any explantation must be reported to your local authorities.

Moreover, in agreement with our local French authority AFSSAPS, Ioltech will continue the follow-up of the patients cohort included in the retrospective study mentioned above, in order to identify significant potential risk factors and/or to adapt the follow-up recommendations for implanted patients.

Ioltech, will continue to keep you updated on this matter. If you have any questions regarding this notice please do not hesitate to contact us directly.

Best Regards,

Pascal Bernard  
Advanced Research Director  
Adverse Event correspondent

Marielle Fournier  
Head of Regulatory Affairs

Contacts:

Pascal Bernard pbernard@ioltech.com  
Tel. +33 5 46 44 08 17

Name of sponsor/Company IOLTECH SA A company of Carl Zeiss Meditec AG	Name of product NEWLIFE (VIVARTE Presbyopic)	Common designation Intra-ocular lens for Presbyopic surgery (anterior chamber)
Title of study	Observational study of tolerance of Newlife implant	
Investigator(s)	Pr Baikoff - Dr Bosc - Dr Simorre - Dr Gauthier - Dr Hachet - Pr Lagoutte / Dr Chastang	
Study centre(s)	6 centres	
Publication (reference)	NA	
Studied period (years) April 2006 – September 2006	Phase of development Post-marketing phase	
Objectives	Following several cases of important loss of endothelial cells after 2-3 year implantation that may lead to corneal damage and explantation of Newlife®. In agreement with local authorities, AFSSAPS, observational study in order to understand the phenomenon and/or risk factors.	
Methodology	Surgeons having implanted the most important numbers of Newlife implant have received a letter explaining the need and objectives of post-operative data analysis. Data have been collected by independent clinical research associates at the centres. Data sheets were completed: <ul style="list-style-type: none"> <li>- pre-operative data</li> <li>- per-operative data</li> <li>- post-operative data at each control visit</li> </ul> When data relative to endothelial cell count were not available at any post-operative visit, patient was not enrolled. An information note has been also provided to the center for the concerned patients, explaining the data collection.	
Number of patients	150 patients (277 eyes) Compliance to eligibility conditions of patients to be implanted with a Newlife has been checked versus surgical protocol recommended by Ioltech.	
Diagnosis and main criteria	See methodology	
Test product	Newlife	
Duration of implantation	Duration of implantation range from 0.26 to 54.47 months (mean = 18.76 months)	
Criteria for evaluation	<ul style="list-style-type: none"> <li>- endothelial cell count, pre-operative and post-operative</li> <li>- variation of endothelial cell count</li> <li>- ACD</li> <li>- Duration of implantation</li> </ul>	
Statistical method	Double data record on Alpha/Open VMS Criteria analysis mainly descriptive. Numeric variables are means, standard deviations, median, 95% confidence interval and extreme values. Qualitative variables are numbers and %.	
Summary/Conclusion	Over the 270 eye data collected, 25 were explanted, among which 6 were explanted due to endothelial cell loss.  Mean duration of implantation observed for the explanted population due to endothelial cell loss is 11 month longer compared to the mean duration of implantation of the global population (29.56 months and 18.76 months respectively).  It seems that a « low » endothelial cell count at pre-operative visit (between 1800 and 2000 c/mm <sup>2</sup> ) is not predictive of pathological cell loss after implantation. <ul style="list-style-type: none"> <li>- none of the explanted case has a pre-operative count below 2000 c/mm<sup>2</sup></li> <li>- none of the global population case having a low pre-operative cell count has been explanted to date for endothelial cell loss and all of them remained stable all over the follow-up.</li> </ul> The potential effect of the diameter of the implant has been also assessed but data collected does not allow to conclude as diameter of the implant	

	<p>and size of irido-corneal angle are both available for only 20.85% of the cases.</p> <p>The depth of anterior chamber seems to be the variable having the most important effect.</p> <p>An ACD below 3.3 mm has been observed in 4 cases of explantations due to endothelial cell loss over 5. e(ACD is not available in the 6th case of explantation due to cell loss).</p> <ul style="list-style-type: none"><li>- Mean ACD in the sub-population explanted due to endothelial cell loss was 3.1 +/- 0.1 mm.</li><li>- Mean ACD in the overall population, excluding the explanted patients is 3.4 +/- 0.3 mm.</li></ul> <p>However, one case of explantation due to endothelial cell loss has an ACD of 3.4 mm that excludes the conclusion that tightness of anterior chamber lead is the single factor leading to dramatic corneal decompensation.</p> <p>These materiovigilance cases may be of multifactor origins. Due to the design of the study, no data is available concerning the way of life of the patients (sport ...).</p> <p>However, extrapolation curves based on the data available so far shows a potential 40-50% cell loss after 60 month implantation.</p> <p>Regular and rigorous survey of corneal endothelium of implanted patients using periodic specular microscopy remains the mean of choice to detect any dramatic corneal decompensation.</p>
--	---