OPHTA-FRANCE S.A.

16 A RUE DE JOUANET  
35700 RENNES

Tél. : 02.99.36.32.33 - Fax. : 02.99.36.82.95

S.A. au capital de 250 020 €

RCS. Rennes 379 437 627 - APE : 4646Z

N° Identification TVA : FR45 379 437 627

Certifié AFAQ - ISO 9001/1996/5756c

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### Description of the problem:

Our contract manufacturer for sterilisation has notified us some deviations in the sterilization process. The sterility of the products is guaranteed but due to these deviations, we cannot exclude that in the concerned batches some lenses could have an increased rate of endotoxins. The affected products could bear an increased risk of immediate postoperative inflammations.

Concerned batches:


- 80 91 52
- 80 91 61
- 80 91 71
- 80 91 82
- 80 92 53

### Corrective action:

Do not to implant the lenses belonging to the batches in the above list and to return it as soon as possible to Xcelens S.A. Those units will be replaced immediately by Xcelens S.A.

### Advise on action to be taken by the user:

Please identify and isolate all lenses belonging to the batches listed above. Your Xcelens S.A. local representative will arrange the exchange of the concerned IOLs.

We remind you that in case of any immediate postoperative inflammation event, the guidelines recommend to proceed to an intensive topical corticosteroids treatment until the inflammation subsides. A brief course of systemic corticosteroids may be beneficial. Frequent follow-up is necessary to monitor IOP and to reassess for signs of bacterial infection.

**Please make sure that this instruction is passed on all those who need to be aware within your organisation or to any customers respectively personnel where the potentially affected device has been transferred.**

In case of need of further information, please contact:

Xcelens S.A.  
Ms. Hélène COCKENPOT  
Tel: +33 4 50 63 50 16  
Fax: +33 4 50 77 66 82  
Email: complaint@xcelens.ch

Xcelens S.A confirms that this notice has been sent to the appropriate National Competent Authorities.

Yours sincerely,

Nikolaus HOFER  
Managing Director  
Xcelens S.A.