Réf. FSN : M18022

Réf FSCA : M18022

Date : 23 octobre 2019

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The User's Manual provided with your equipment does not include all the Maintenance information for the device.

Please insert the following pages at the beginning of chapter Maintenance of your User's Manual.

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## User Manual Addendum

The expected service life for Apollo angulation chain is 8 years.

Contact your service organization in order to plan the intervention once the unit is close to reach the next replacement date given in the table below.

At replacement, FSE will fill the column “Next replacement date” adding 8 years to the replacement date in both User and Service Manual

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The Service Manual provided with your equipment does not include all the Maintenance information for the device.

Please insert the following pages at the beginning of chapter “Maintenance and Cleaning” of your Service Manual.

All other specifications are unchanged with respect to the ones provided in your Manual.

## Service Manual Addendum

### 9.2 Parts to be replaced periodically

The expected lifetime of Apollo angulation chain is 8 years. After the operation time of 8 years, the following part must be replaced:

Code	Description	Replace after	Installation date
6690702200	Angulation chain	Every 8 years	

At replacement, FSE have to fill the column “Next replacement date” adding 8 years to the replacement date in both User and Service Manual

Angulation chain replacement date (dd – mm – yyyy)	Reason of the replacement  corrective maintenance / planned maintenance	next replacement date mm/yyyy	Executed by Organization (FSE Name)	Date / Signature

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