AVIS DE SECURITE URGENT

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URGENT FIELD SAFETY NOTICE

Commercial name of the affected product: T4P

REF : 1000587

FSCA-identifier: xx

Type of action: retrofit (battery exchange, and software upgrade)

Date: 26 February 2020

Attention:

Technical managers and technicians of home health care providers using T4P device.

DETAILS ON AFFECTED PRODUCT

TYPE OF DEVICE: Air pressure transducer

MODEL NAME: T4P

REFERENCE (REF): 1000587

BATCH (LOT)	DATE OF MANUFACTURE:
18L7856	fabriqué semaine 24 2018 ;
18L7857	fabriqué semaine 26 2018 ;
18L7860	fabriqué semaine 28 2018 ;
18L7861	fabriqué semaine 36 2018.

Dear Madam, dear sir,

SRETT implements a field safety corrective action. The purpose of this letter is to inform you of the associated risks and the actions you will need to take. I kindly request that you confirm receipt of this letter using the form in Appendix I.

The T4P devices, concerned by this notice are **from the following batches:**

18L7856 manufactured week 24 2018; **18L7857** manufactured week 26 2018;

18L7860 manufactured week 28 2018; **18L7861** manufactured week 36 2018.

Thanks to the traceability of the devices, we are able to provide a **list of the devices** concerned in your possession, with their current state. We will send you this list by email. You will find in appendix II how to read batch number and serial number of you T4P devices

DESCRIPTION OF THE PROBLEM

We found that the batteries in the T4P devices concerned are depleting faster than specified, and faster than indicated on the remote monitoring interface. Under these conditions, the devices may stop working before the user is instructed to replace the batteries.

The user manual indicates a three-year replaceable battery life. Our analysis indicates that the actual battery life is about a year and a half, with a variability of several months.

Before the battery runs out, devices are operational and the data transmitted are complete and valid.

However, once the battery is depleted, the device ceases to operate, to transmit data, and its indicator lights no longer function. The absence of transmitted data is indicated by a "transmission failure" indicator in the telemonitoring interface.

We have also found that an estimated 30% of devices that have stopped working after their battery has run down, do not restart without exceptional operation, after battery replacement.

Impact of patient's safety

There is no risk of direct impact on patient's safety from using one of the devices according to the instructions in the user manual even if the battery is exhausted. In particular, there are no consequences on the application of the treatment being carried out by the continuous positive pressure machine.

Origin of the fault and correction

The lack of autonomy of the devices is due to a malfunction of the embedded software. This malfunction creates an electrical over-consumption. This over-consumption has no other impact than accelerating battery depletion.

We have identified this malfunction and are working on the development of a new version of the software. The software update in the devices will be done remotely through their GPRS connectivity, while the devices are in operation.

Nevertheless, it is necessary to change the battery in the devices, even if the device continues to work, as device battery has been discharged prematurely.

ADVISE ON ACTION TO BE TAKEN BY THE USER

We propose the measures described below. We would like to accompany you in their implementation and are at your disposal for this purpose. We kindly ask the technical teams to carry out the operations described below, as soon as possible.

A. For the non-deployed devices from the batches concerned

Isolate the devices, and quarantine them while waiting for the new software version to be released. We will let you know when the new version is available, so that you can change the battery.

B. For devices deployed in patients from the batches concerned

1. Change the battery of the devices concerned according to the procedure described in the attached instruction sheet "1002705.AA T4P - Battery replacement procedure"
2. Check their correct operation after changing the battery according to the procedure described in the attached instruction sheet "1002706.AA T4P - Activation Procedure".
3. For devices that were found to be non-functional in Step 2, return them to SRETT so that we can perform the technical operation needed to restart them. We will then, return them

to you. We are attaching to this letter, a form, « 1001180.AB Fiche_retour_SAV » , to document product returns.

Timeline

By changing the battery quickly, you will help us to avoid an interruption in the operation of the devices.

We suggest that we hold without delay a meeting to organize together the operations described above and to plan our support.

TRANSMISSION OF THIS FIELD SAFETY NOTICE

This notice needs to be passed on all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred.

Please transfer this notice to other organizations on which this action has an impact.

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This Safety Notice has been communicated to the following National Competent Authority: ANSM, National Agency for the Safety of Medicines and Health Products (France)

Please accept our apologies for the inconvenience caused by this corrective action. This voluntary corrective action reaffirms our commitment to develop, manufacture and market devices of the highest quality for the remote monitoring of patients in home care. We are grateful for your cooperation in this process.

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ANNEXE I / APPENDIX I
CONFIRMATION DE RECEPTION / ACKNOWLEDGEMENT OF RECEIPT

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We kindly request the designated person in your organization to return this acknowledgement of receipt of this field safety notice.

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