# MassCheck<sup>®</sup> Neuroleptics 1/*EXTENDED* Plasma Control Level II

Order no. 0212/XT

## Instructions

These Chromsystems plasma controls are designed to monitor the accuracy and precision of analytical procedures for the quantitative determination of *neuroleptics* in serum/plasma.

The lyophilised control is based on human plasma and should be used routinely according to the instructions of the assay procedure being used. While following these instructions the control is to be handled in the same manner as a patient specimen.

## Reconstitution

Carefully remove the rubber stopper from the vial. Add exactly 1.0 ml of distilled water into the vial. Replace the stopper and allow to sit for 10 to 15 minutes at room temperature. Swirl the vial to dissolve the contents until homogeneity.

## Storage conditions

Stored below -18 °C, unopened controls are stable until the expiry date which is stated on the vial labels. Light protected and tightly capped, the reconstituted controls are stable for 14 days at +2 to +8 °C except Olanzapine\*\* and Desmethylolanzapine\*\* which should be aliquoted immediately after reconstitution and stored below -18 °C. If the controls cannot be analysed within this period, they should be aliquoted and stored frozen below -18 °C (up to 3 months).

## Value Assignment

The mean values and acceptable ranges have been determined from a representative sampling of this control lot. These values are specific for this control lot. The individual laboratory should establish their own mean values and acceptable ranges based on instrument, reagent and laboratory variations. Due to procedural and reagent modifications the mean values and acceptable ranges may alter slightly.

## Caution

This product is manufactured from pooled human plasma which has been tested by the manufacturer and found negative for infections by the human immunodeficiency virus (HIV), the hepatitis B virus (HBV), the hepatitis C virus (HCV) and the bacterium *Treponema pallidum*. Nevertheless, a potential risk of infection cannot be entirely excluded. Consider all products containing human source material as potentially infectious and exercise the same care in the handling of this product as in the handling of potentially infectious patient samples.

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Max. variance flacon à flacon: 1 % (poids)

## Concentrations

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**Concentrations/Concentrations**  
**MassCheck® Neuroleptics 1/EXTENDED**  
**Plasma Control Level II**

**Order no/N° de commande: 0212/XT**

**Lot no: 2219**

**Exp. Date/Date d'expiration: 2022-05-01**

<b>Substance</b>	<b>Unit</b>	<b>Target value*</b>	<b>Range</b>	
<b>Substance</b>	<b>Unité</b>	<b>Valeur cible*</b>	<b>Intervalle</b>	
Aripiprazole	µg/l	304	243	- 365
Dehydroaripiprazole	µg/l	122	97,9	- 147
Clozapine	µg/l	469	375	- 562
Desmethylclozapine	µg/l	204	163	- 245
Haloperidol	µg/l	18,0	14,4	- 21,6
Olanzapine**	µg/l	68,0	54,4	- 81,6
N-Desmethylolanzapine**	µg/l	73,1	58,5	- 87,8
Quetiapine	µg/l	263	211	- 316
Norquetiapine	µg/l	96,3	77,0	- 116
Risperidone	µg/l	15,2	12,2	- 18,2
9-OH-Risperidone	µg/l	64,1	51,3	- 76,9

**\*Method of Analysis/Méthode d'analyse:**

Chromsystems reagent kit **MassTox®** TDM BASIC Kit Series A, order no. / N° der commande 92111 and/et **MassTox®** TDM PARAMETER-Set Neuroleptics 1/EXTENDED in serum/plasma, order no. / N° der commande 92912/XT

**Traceability/Traçabilité**

The values of the control level of Dehydroaripiprazole, Clozapine, N-Desmethylclozapine, Haloperidol, Olanzapine, Quetiapine, Norquetiapine, Risperidone and 9-OH-Risperidone are traceable to certified reference materials and primary standards.

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**\*\*Note:**

These analytes are not stable at +2 to +8 °C and at room temperature. The controls should be aliquoted immediately after reconstitution and stored below -18 °C.

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