Boulogne Billancourt,

Le 15 décembre 2010

### INFORMATION PRODUIT IMPORTANTE

Sondes endocavitaires St. Jude Medical Riata et Riata ST Silicone  
Modèles 1560, 1561, 1562, 1570, 1571, 1572, 1580, 1581, 1582, 1590, 1591, 1592,  
7000, 7001, 7002, 7010, 7011, 7040, 7041, 7042

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Kathleen M. Chester  
Sr. Vice President, Regulatory Affairs and Quality Assurance

## REFERENCES

1. Mehta D, Nayak H, Singson M, Chao S, Jorge E, Camunas L, Gomes J. Late Complications in Patients with Pectoral Defibrillator Implants with Transvenous Defibrillator Lead Systems: High Incidence of Insulation Breakdown. *PACE*, 1998; 21:1893-1900.
2. De Lurgio D, Sathavorn C, Mera F, Leon A, Walter P, Langberg J. Incidence and Implications of Abrasion of Implantable Cardioverter-Defibrillator Leads. *American Journal of Cardiology*, 1997, 79:1409-1411.
3. Pavia S, Saliba W, Wilkoff B. Lead System Dysfunction, Diagnosis, and Therapy. In: Pacifico A. ed-in-chief. *Implantable Defibrillator Therapy: A Clinical Guide*. Norwell, Ma. Kluwer Academic Publishers 2002:259-278.
4. Kleemann T, Becker T, Doenges K, Vater M, Senges J, Schneider S, Saggau W, Weisse U, Seidl K. Annual Rate of Transvenous Defibrillation Lead Defects in Implantable Cardioverter-Defibrillators over a Period > 10 years. *Circulation*. 2007;115:2474-2480.
5. Blommaert D, Louagie Y, Eucher P, Collet B, Mancini I, Govarets G, Leclercq C, De Roy L. Are Cardioverter-Defibrillator Leads Reliable Over time? *NASPE Abstracts*, 2003, 26:1123.
6. Degeratu F, Khalighi K, Peters R, Shorofsky S, Gold M. Sensing Lead Failure in Implantable Defibrillators: A Comparison of Two Commonly Used Leads. *J Cardiovasc Electrophysio*, 2000, 11:21-24.
7. Valk S, Luijten R, Jordaens L. Insulation Damage in a Shock Wire: An Unexpected Fluoroscopic Image. *PACE*. 2010; 1-3.
8. Richards M, Warren C, Anderson M. Late failure of a single-coil Transvenous implantable cardioverter-defibrillator lead associated with conductor separation. *Europace Advanced Access* published March 27, 2010.
9. Jalal Z, Derval N, Ploux S, Bordachar P. Unusual failure of a multilumen, small-diameter implantable cardioverter-defibrillator lead. *Heart Rhythm*, 2009, Articles in Press.
10. Duray G, Israel C, Schmitt J, Hohnloser S. Implantable cardioverter-defibrillator lead disintegration at that level of the tricuspid valve. *Heart Rhythm Society*, 2008.
11. HRS/EHRA Expert Consensus on Monitoring Cardiovascular Implantable Electronic Devices (CIED) dated April 2008
12. Wilkoff et al., Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management, *Heart Rhythm*, Vol 6, No 7, July 2009.
13. Maisel et al., Recommendations from the Heart Rhythm Society Task Force on Lead Performance Policies and Guidelines, *Heart Rhythm*, Vol 6, No 6, June 2009.