**-RAPPEL DE LOTS-
URGO Traitement Taches Brunes
Code 551849**

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Cher client, cher confrère,

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Chenôve, July 10th 2017,

Dear client, Dear colleague,

Following one event in which a user got a superficial burn after use of the Urgo Brown Spot Treatment Pen, preliminary conclusions regarding this event have led us to believe that the product's packaging and instruction for use need to be improved. All necessary information for safe use of the product is currently provided in the instruction for the user and outer box, but our risk assessment has shown that the information is not appealing to read to users, which could lead to skin impact due to use error.

Taking into account the possible risk of non full observance of the instructions for use, and thus the risk of a superficial burn with new users, Laboratories Urgo decided by precaution, in accordance with The French Authority ANSM and with the manufacturer Medical Brands, to recall the product Urgo Traitement Taches Brunes (CIP : 3401360210110).

Concerned batches:

610320 / Expiry date : 2018/09 – Mfg date : 2016/09

610321 / Expiry date : 2018/09 - Mfg date : 2016/09

611300 / Expiry date : 2018/09 - Mfg date : 2016/09

This recall is limited to pharmacists and distributors, since the risk assessment has shown that the product can still be used safely if used according to the instructions for use.

The medical device is developed, manufactured and packaged by the company Medical Brands. The distribution in France is managed by Laboratoires URGO, under the brand name URGO Traitement Taches Brunes.

Laboratoires URGO will contact you within 48 hours for the conditions of return and reimbursement.

We can confirm that this product recall procedure has been communicated to the ANSM and we are sorry for the inconvenience.

Best Regards,

Marie-Catherine Sarciaux
Qualified Person

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Medical Brands B.V.