
Avis de sécurité urgent sur le terrain

Nom commercial / modèle: Tampons de gaze avec fil à rayons X

Identifiant FSCA : FSCA2020001

Type d'action : Rappel

2020-08-18

Attention : [[Medline Int'l France SAS](#)]

Cher Monsieur ou Madame,

Cette lettre vous fournira les informations suivantes :

Détails sur les dispositifs concernés : Tampons de gaze avec fil à rayons X

Veillez vous référer à l'annexe- [2&3-Number of Medical Device Concern to Recall in Healthcare update](#)

Description du problème :

Sur la base des commentaires de notre client, nous avons appris que les fils de contraste pour rayons X sous tension dans les compresses de gaze avec un fil de contraste pour rayons X peuvent se casser et se briser à des points individuels. Veuillez vous référer l'annexe- [2&3-Number of Medical Device Concern to Recall in Healthcare update](#) pour les informations des dispositifs concernés.

Nous informons le Medline Int'l France SAS pour isoler les dispositifs concernés dans l'entrepôt.

En général, de petits morceaux de fil peuvent se produire lorsque le fil de rayons X se casse ou se brise. Cela pourrait, par exemple, conduire à une inflammation et / ou à la formation de granulomes en restant dans le corps.

Conseiller sur les mesures à prendre par l'administrateur de l'hôpital :

1. Veuillez transmettre cet avis à tous ceux qui doivent en être informés dans votre organisation ou n'importe où le ou les dispositif(s) potentiellement concerné(s) ont été transférés.
2. Veuillez arrêter de les utiliser, nous rappellerons les produits concernés qui ont été vendus mais non utilisés.

Conseiller sur les mesures à prendre par le distributeur :

1. Veuillez transmettre cet avis à tous ceux qui doivent en être informés dans votre organisation ou n'importe où le ou les dispositif(s) potentiellement concerné(s) ont été délivrés...
2. Isolez immédiatement les dispositifs concernés situés dans les entrepôts.
3. Veuillez arrêter de les vendre, nous rappellerons les produits concernés qui ont été vendus mais non utilisés.

Transmission de cet avis de sécurité sur le terrain :

Cet avis doit être transmis à tous ceux qui doivent en être informés dans votre organisation ou n'importe où le ou les dispositif(s) potentiellement concerné(s) ont été transférés.

Veuillez rester au courant de cet avis et des mesures qui en découlent pendant une période appropriée pour garantir l'efficacité de l'action corrective.

Nous vous serions reconnaissants de bien vouloir confirmer la réception de cette lettre. Veuillez remplir le formulaire d'accusé de réception ci-dessous et le renvoyer par e-mail ou fax.

Contacts de référence :

Nous vous prions de nous excuser pour cette inconvenance. Si vous avez des questions, veuillez contacter votre directeur général local de la qualité - Vince Tian

Organisation : Allmed Medical Products Co.,Ltd

Tél : 0086- (0) 717- 4215906

Fax : 0086- (0) 717-4215989

E-mail : vince@allmed.cn

Cet avis a été notifié à l'agence de réglementation appropriée.

(Paragraphe de conclusion)

Signature :

.....

Directeur général, Centre de qualité

No.18 Qixing Road, Majiadian Town, 443200 Zhijiang City, Hubei Province, Chine

Tél : 0086-(0)717- 4215906

Fax : 0086-(0)717-4215989

E-mail : vince@allmed.cn

Formulaire d'accusé de réception

=====
Confirmation de réception de l'avis de sécurité sur le terrain
Produits concernés : veuillez vous référer à l'annexe- 2&3-Number of Medical Device
Concern to Recall in Healthcare update
FSCA : FSCA2020001
Type de FSCA : Rappel
=====

Veillez remplir ce formulaire et renvoyer immédiatement cette confirmation par e-mail ou fax.

Fax : 0086-(0)717-4215989

E-mail : vince@allmed.cn

Nom : Vince Tian

Tél. No. : 0086-(0)717- 4215906

Adresse e-mail : vince@allmed.cn

Date et signature : _____

Adresse de l'organisation :

 No.18 Qixing Road, Majiadian Town, 443200Zhijiang City, Hubei Province, Chine

En tant que distributeur des produits, ainsi que les produits sont les marques de notre client Medline, les informations suivantes seront collectées par Medline.

1. Liste des... distribués avec des équipements

Pays / Région	Équipement S / N	Nom commercial / Modèle	Distributeur / Hôpital	Contact	Adresse	Téléphone	E-mail

2. Liste des....individuels distribués

Pays / Région	Nom commercial / Modèle	Quantité concernée	Distributeur / Hôpital	Contact	Adresse	Téléphone	E-mail



ALLMED MEDICAL PRODUCTS CO.,LTD
No.18 Qi Xing Road, Majiadian Town, 443200
Zhijiang City, Hubei Province, P.R.C
Tel: 86 717 4211111 Fax: 86 717 4225499

Number of Medical Device Concern to Recall in Healthcare

Date: August 21,2020

Dear Sir/Madam from BfArM:

Allmed Medical has received NCAR 11893/20 from BfArM regarding FSCA2020001 on Aug14 2020,for response the email in a timely manner, we prepare detailed product list including all requirements that NCAR need as follows. Allmed Medical send this list on Aug19, 2020 and update it on Aug21, 2020 for your reference.

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
1	Laboratoire Tetra Medical	77492	Gauze Swabs with X-ray Thread	01039K2	PKG	4000	40	1.223	2019/9/8	A523-19080380	774920003	France	Compresse Gaze
2	Laboratoire Tetra Medical	77353	Gauze Swabs with X-ray Thread	01039K2	PKG	4000	40	1.223	2019/8/29	A523-19080312	773530003	France	Compresse Gaze
3	Laboratoire Tetra Medical	77683	Gauze Swabs with X-ray Thread	01039K2	PKG	4000	40	1.223	2019/11/14	A523-19110012	776830004	France	Compresse Gaze
4	Laboratoire Tetra Medical	77685/77972/77823/ 77967/77970/	Gauze Swabs with X-ray Thread	01039K2	PKG	2000	20	0.611	2020/3/8	A523-20020162	779700004	France	Compresse Gaze
5	Laboratoire Tetra Medical	78650	Gauze Swabs with X-ray Thread	22113K2	PKG	5400	36	1.197	2020/6/10	A523-20060138	786500005	France	Compresse Gaze
6	Laboratoire Tetra Medical	78650	Gauze Swabs with X-ray Thread	22120K2	PKG	800	10	0.342	2020/6/10	A523-20060138	786500006	France	Compresse Gaze
7	Laboratoire Tetra Medical	76983	Gauze Swabs with X-ray Thread	01173T2	PKG	7200	36	1.508	2019/7/2	A523-19070035	769830006	France	Compresse Gaze
8	Laboratoire Tetra Medical	77177	Gauze Swabs with X-ray Thread	01173T2	PKG	14400	72	3.016	2019/7/31	A523-19070435	771770009	France	Compresse Gaze
9	Laboratoire Tetra Medical	77360	Gauze Swabs with X-ray Thread	01173T2	PKG	34000	170	7.121	2019/8/8	A523-19080065	773600004	France	Compresse Gaze
10	Laboratoire Tetra Medical	77681	Gauze Swabs with X-ray Thread	01173T2	PKG	64800	324	13.57	2019/11/4	A523-19100308	776810007	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
11	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	01173T2	PKG	35000	175	7.331	2020/1/6	A523-20010019	778210004	France	Compresse Gaze
12	Laboratoire Tetra Medical	76985/76981	Gauze Swabs with X-ray Thread	01172T2	PKG	25760	322	13.26	2019/7/2	A523-19070049	769850002	France	Compresse Gaze
13	Laboratoire Tetra Medical	76983	Gauze Swabs with X-ray Thread	01172T2	PKG	57840	723	29.77	2019/7/2	A523-19070035	769830005	France	Compresse Gaze
14	Laboratoire Tetra Medical	76976/76798/77158	Gauze Swabs with X-ray Thread	01172T2	PKG	36800	460	18.94	2019/7/8	A523-19070037	769760003	France	Compresse Gaze
15	Laboratoire Tetra Medical	77174/76983	Gauze Swabs with X-ray Thread	01172T2	PKG	560	7	0.288	2019/7/12	A523-19070215	769830005	France	Compresse Gaze
16	Laboratoire Tetra Medical	77177	Gauze Swabs with X-ray Thread	01172T2	PKG	52800	660	27.18	2019/7/31	A523-19070435	771770008	France	Compresse Gaze
17	Laboratoire Tetra Medical	77168/77169	Gauze Swabs with X-ray Thread	01172T2	PKG	48000	600	24.71	2019/7/31	A523-19070412	771680003	France	Compresse Gaze
18	Laboratoire Tetra Medical	77494	Gauze Swabs with X-ray Thread	01172T2	PKG	32000	400	16.47	2019/9/8	A523-19090051	774940002	France	Compresse Gaze
19	Laboratoire Tetra Medical	77495/77532	Gauze Swabs with X-ray Thread	01172T2	PKG	38000	475	19.56	2019/9/16	A523-19090116	774950001	France	Compresse Gaze
20	Laboratoire Tetra Medical	77604/77398	Gauze Swabs with X-ray Thread	01172T2	PKG	37440	468	19.27	2019/10/14	A523-19100025	776040003	France	Compresse Gaze
21	Laboratoire Tetra Medical	77682/77717	Gauze Swabs with X-ray Thread	01172T2	PKG	72000	900	37.06	2019/11/14	A523-19110013	776820005	France	Compresse Gaze
22	Laboratoire Tetra Medical	77814	Gauze Swabs with X-ray Thread	01172T2	PKG	7200	90	3.706	2020/1/1	A523-20010014	778140005	France	Compresse Gaze
23	Laboratoire Tetra Medical	77822	Gauze Swabs with X-ray Thread	01172T2	PKG	67840	848	34.92	2020/1/1	A523-20010013	778220003	France	Compresse Gaze
24	Laboratoire Tetra Medical	78316/78318/78199	Gauze Swabs with X-ray Thread	01172T2	PKG	70000	875	36.03	2020/4/24	A523-20040485	781990002	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
25	Laboratoire Tetra Medical	78643/78719	Gauze Swabs with X-ray Thread	01172T2	PKG	72000	900	37.06	2020/5/10	A523-20050161	786430004	France	Compresse Gaze
26	Laboratoire Tetra Medical	76983	Gauze Swabs with X-ray Thread	20609T2	PKG	6000	240	5.431	2019/7/2	A523-19070035	769830008	France	Compresse Gaze
27	Laboratoire Tetra Medical	77176	Gauze Swabs with X-ray Thread	20609T2	PKG	15600	624	14.12	2019/7/19	A523-19070248	771760009	France	Compresse Gaze
28	Laboratoire Tetra Medical	77601	Gauze Swabs with X-ray Thread	20609T2	PKG	15600	624	14.12	2019/10/5	A523-19090299	776010010	France	Compresse Gaze
29	Laboratoire Tetra Medical	77681	Gauze Swabs with X-ray Thread	20609T2	PKG	1200	48	1.086	2019/11/4	A523-19100308	776810008	France	Compresse Gaze
30	Laboratoire Tetra Medical	77820/77687	Gauze Swabs with X-ray Thread	20609T2	PKG	7200	288	6.518	2020/1/1	A523-19120451	778200008	France	Compresse Gaze
31	Laboratoire Tetra Medical	78316/78180/78181/ 78196/78197/	Gauze Swabs with X-ray Thread	20609T2	PKG	9000	360	8.147	2020/4/3	A523-20030481	783160010	France	Compresse Gaze
32	Laboratoire Tetra Medical	78180/78197/78199/ 78309/78313/	Gauze Swabs with X-ray Thread	20609T2	PKG	17000	680	15.39	2020/4/8	A523-20040197	781990004	France	Compresse Gaze
33	Laboratoire Tetra Medical	78722	Gauze Swabs with X-ray Thread	20609T2	PKG	9600	384	8.69	2020/6/6	A523-20050609	787220003	France	Compresse Gaze
34	Laboratoire Tetra Medical	76985/76981	Gauze Swabs with X-ray Thread	22132T2	PKG	8400	168	3.156	2019/7/2	A523-19070049	769850003	France	Compresse Gaze
35	Laboratoire Tetra Medical	77176	Gauze Swabs with X-ray Thread	22132T2	PKG	2800	56	1.052	2019/7/19	A523-19070248	771760010	France	Compresse Gaze
36	Laboratoire Tetra Medical	77682/77717	Gauze Swabs with X-ray Thread	22132T2	PKG	2800	56	1.052	2019/11/14	A523-19110013	776820006	France	Compresse Gaze
37	Laboratoire Tetra Medical	77820/77687	Gauze Swabs with X-ray Thread	22132T2	PKG	2800	56	1.052	2020/1/1	A523-19120451	778200009	France	Compresse Gaze
38	Laboratoire Tetra Medical	78628	Gauze Swabs with X-ray Thread	22132T2	PKG	2800	56	1.052	2020/5/10	A523-20050159	786280008	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
39	Laboratoire Tetra Medical	76983	Gauze Swabs with X-ray Thread	20606T2	PKG	16800	560	12.2	2019/7/2	A523-19070035	769830007	France	Compresse Gaze
40	Laboratoire Tetra Medical	77175	Gauze Swabs with X-ray Thread	20606T2	PKG	28800	960	20.91	2019/7/17	A523-19070214	771750007	France	Compresse Gaze
41	Laboratoire Tetra Medical	77176	Gauze Swabs with X-ray Thread	20606T2	PKG	21600	720	15.68	2019/7/19	A523-19070248	771760008	France	Compresse Gaze
42	Laboratoire Tetra Medical	77597	Gauze Swabs with X-ray Thread	20606T2	PKG	8400	280	6.098	2019/9/30	A523-19090305	775970007	France	Compresse Gaze
43	Laboratoire Tetra Medical	77684	Gauze Swabs with X-ray Thread	20606T2	PKG	20160	672	14.63	2019/12/6	A523-19120093	776840006	France	Compresse Gaze
44	Laboratoire Tetra Medical	77820/77687	Gauze Swabs with X-ray Thread	20606T2	PKG	10080	336	7.317	2020/1/1	A523-19120451	778200007	France	Compresse Gaze
45	Laboratoire Tetra Medical	78718/78710	Gauze Swabs with X-ray Thread	20606T2	PKG	16800	560	12.2	2020/6/10	A523-20060139	787100007	France	Compresse Gaze
46	Laboratoire Tetra Medical	78644/78783	Gauze Swabs with X-ray Thread	20606T2	PKG	25200	840	18.29	2020/5/10	A523-20050160	786440007	France	Compresse Gaze
47	Laboratoire Tetra Medical	76983	Gauze Swabs with X-ray Thread	20612T2	PKG	1560	78	1.433	2019/7/2	A523-19070035	769830009	France	Compresse Gaze
48	Laboratoire Tetra Medical	77174/76983	Gauze Swabs with X-ray Thread	20612T2	PKG	1560	78	1.433	2019/7/12	A523-19070215	771740011	France	Compresse Gaze
49	Laboratoire Tetra Medical	77603	Gauze Swabs with X-ray Thread	20612T2	PKG	1560	78	1.433	2019/10/14	A523-19100024	776030015	France	Compresse Gaze
50	Laboratoire Tetra Medical	77681	Gauze Swabs with X-ray Thread	20612T2	PKG	1560	78	1.433	2019/11/4	A523-19100308	776810009	France	Compresse Gaze
51	Laboratoire Tetra Medical	78189/78190/78192/ 78193/78183/	Gauze Swabs with X-ray Thread	20612T2	PKG	4680	234	4.3	2020/3/18	A523-20030173	781920008	France	Compresse Gaze
52	Laboratoire Tetra Medical	78720	Gauze Swabs with X-ray Thread	20612T2	PKG	4500	225	4.134	2020/5/23	A523-20050348	787200005	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
53	Laboratoire Tetra Medical	76976/76798/77158	Gauze Swabs with X-ray Thread	20608T2	PKG	3060	102	2.794	2019/7/8	A523-19070037	769760004	France	Compresse Gaze
54	Laboratoire Tetra Medical	77180	Gauze Swabs with X-ray Thread	20608T2	PKG	3060	102	2.794	2019/8/1	A523-19070453	771800003	France	Compresse Gaze
55	Laboratoire Tetra Medical	77361	Gauze Swabs with X-ray Thread	20608T2	PKG	1530	51	1.397	2019/8/15	A523-19080128	773610010	France	Compresse Gaze
56	Laboratoire Tetra Medical	77493/77354	Gauze Swabs with X-ray Thread	20608T2	PKG	1530	51	1.397	2019/9/8	A523-19090045	774930007	France	Compresse Gaze
57	Laboratoire Tetra Medical	77683	Gauze Swabs with X-ray Thread	20608T2	PKG	4590	153	4.191	2019/11/14	A523-19110012	776830012	France	Compresse Gaze
58	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	20608T2	PKG	1530	51	1.397	2020/1/6	A523-20010019	778210009	France	Compresse Gaze
59	Laboratoire Tetra Medical	78726	Gauze Swabs with X-ray Thread	20608T2	PKG	3060	102	2.794	2020/6/16	A523-20060284	787230008	France	Compresse Gaze
60	Laboratoire Tetra Medical	76976/76798/77158	Gauze Swabs with X-ray Thread	20701T2	PKG	9000	180	8.051	2019/7/8	A523-19070037	769760002	France	Compresse Gaze
61	Laboratoire Tetra Medical	77169	Gauze Swabs with X-ray Thread	20701T2	PKG	7500	150	6.71	2019/7/28	A523-19070345	771690002	France	Compresse Gaze
62	Laboratoire Tetra Medical	77494	Gauze Swabs with X-ray Thread	20701T2	PKG	13500	270	12.08	2019/9/8	A523-19090051	774940001	France	Compresse Gaze
63	Laboratoire Tetra Medical	77603	Gauze Swabs with X-ray Thread	20701T2	PKG	6000	120	5.368	2019/10/14	A523-19100024	776030010	France	Compresse Gaze
64	Laboratoire Tetra Medical	77680	Gauze Swabs with X-ray Thread	20701T2	PKG	4500	90	4.026	2019/10/24	A523-19100242	776800007	France	Compresse Gaze
65	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	20701T2	PKG	4500	90	4.026	2020/1/6	A523-20010019	778210001	France	Compresse Gaze
66	Laboratoire Tetra Medical	78316/78180/78181/ 78196/78197/	Gauze Swabs with X-ray Thread	20701T2	PKG	13000	260	11.63	2020/4/3	A523-20030481	781970003	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
67	Laboratoire Tetra Medical	76976/76798/77158	Gauze Swabs with X-ray Thread	22134T2	PKG	5760	72	2.658	2019/7/8	A523-19070037	769760005	France	Compresse Gaze
68	Laboratoire Tetra Medical	77361	Gauze Swabs with X-ray Thread	22134T2	PKG	5760	72	2.658	2019/8/15	A523-19080128	773610011	France	Compresse Gaze
69	Laboratoire Tetra Medical	77494	Gauze Swabs with X-ray Thread	22134T2	PKG	2880	36	1.329	2019/9/8	A523-19090051	774940004	France	Compresse Gaze
70	Laboratoire Tetra Medical	77680	Gauze Swabs with X-ray Thread	22134T2	PKG	2880	36	1.329	2019/10/24	A523-19100242	776800013	France	Compresse Gaze
71	Laboratoire Tetra Medical	77969/77970	Gauze Swabs with X-ray Thread	22134T2	PKG	2880	36	1.329	2020/2/22	A523-20020064	779700009	France	Compresse Gaze
72	Laboratoire Tetra Medical	77174/76983	Gauze Swabs with X-ray Thread	01167T2	PKG	10800	36	1.42	2019/7/12	A523-19070215	771740009	France	Compresse Gaze
73	Laboratoire Tetra Medical	77358/77351	Gauze Swabs with X-ray Thread	01167T2	PKG	10800	36	1.42	2019/8/4	A523-19080007	773580014	France	Compresse Gaze
74	Laboratoire Tetra Medical	77685	Gauze Swabs with X-ray Thread	01167T2	PKG	10800	36	1.42	2019/12/13	A523-19120172	776850010	France	Compresse Gaze
75	Laboratoire Tetra Medical	78949	Gauze Swabs with X-ray Thread	01167T2	PKG	10800	36	1.42	2020/7/14	A523-20070183	789490008	France	Compresse Gaze
76	Laboratoire Tetra Medical	78644/78783	Gauze Swabs with X-ray Thread	01167T2	PKG	10800	36	1.42	2020/5/10	A523-20050160	786440005	France	Compresse Gaze
77	Laboratoire Tetra Medical	77174/76983	Gauze Swabs with X-ray Thread	01191T2	PKG	18000	90	4.664	2019/7/12	A523-19070215	771740010	France	Compresse Gaze
78	Laboratoire Tetra Medical	77358/77351	Gauze Swabs with X-ray Thread	01191T2	PKG	6000	30	1.555	2019/8/4	A523-19080007	773580015	France	Compresse Gaze
79	Laboratoire Tetra Medical	77490	Gauze Swabs with X-ray Thread	01191T2	PKG	12000	60	3.11	2019/8/29	A523-19080326	774900009	France	Compresse Gaze
80	Laboratoire Tetra Medical	77601	Gauze Swabs with X-ray Thread	01191T2	PKG	12000	60	3.11	2019/10/5	A523-19090299	776010009	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
81	Laboratoire Tetra Medical	77685	Gauze Swabs with X-ray Thread	01191T2	PKG	12000	60	3.11	2019/12/13	A523-19120172	776850011	France	Compresse Gaze
82	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	01191T2	PKG	6000	30	1.555	2020/1/6	A523-20010019	778210006	France	Compresse Gaze
83	Laboratoire Tetra Medical	78723	Gauze Swabs with X-ray Thread	01191T2	PKG	12000	60	3.11	2020/6/10	A523-20060136	787230014	France	Compresse Gaze
84	Laboratoire Tetra Medical	77175	Gauze Swabs with X-ray Thread	01081T2	PKG	10800	360	8.026	2019/7/17	A523-19070214	771750006	France	Compresse Gaze
85	Laboratoire Tetra Medical	77182	Gauze Swabs with X-ray Thread	01081T2	PKG	10800	360	8.026	2019/7/31	A523-19070439	771820005	France	Compresse Gaze
86	Laboratoire Tetra Medical	77680	Gauze Swabs with X-ray Thread	01081T2	PKG	3600	120	2.675	2019/10/24	A523-19100242	776800010	France	Compresse Gaze
87	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	01081T2	PKG	3600	120	2.675	2020/1/6	A523-20010019	778210002	France	Compresse Gaze
88	Laboratoire Tetra Medical	78642	Gauze Swabs with X-ray Thread	01081T2	PKG	1800	60	1.338	2020/4/30	A523-20040581	786420006	France	Compresse Gaze
89	Laboratoire Tetra Medical	77178	Gauze Swabs with X-ray Thread	01174T2	PKG	2E+05	1890	61.73	2019/7/19	A523-19070261	771780002	France	Compresse Gaze
90	Laboratoire Tetra Medical	77495/77532	Gauze Swabs with X-ray Thread	01174T2	PKG	1E+05	1440	47.03	2019/9/16	A523-19090116	774950002	France	Compresse Gaze
91	Laboratoire Tetra Medical	77397/77496	Gauze Swabs with X-ray Thread	01174T2	PKG	36000	450	14.7	2019/9/28	A523-19090263	774960003	France	Compresse Gaze
92	Laboratoire Tetra Medical	77604/77398	Gauze Swabs with X-ray Thread	01174T2	PKG	14480	181	0	2019/10/14	A523-19100025	776040004	France	Compresse Gaze
93	Laboratoire Tetra Medical	77605	Gauze Swabs with X-ray Thread	01174T2	PKG	56000	700	22.86	2019/10/27	A523-19100208	776050004	France	Compresse Gaze
94	Laboratoire Tetra Medical	77683	Gauze Swabs with X-ray Thread	01174T2	PKG	18400	230	7.512	2019/11/14	A523-19110012	776830008	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
95	Laboratoire Tetra Medical	77678	Gauze Swabs with X-ray Thread	01174T2	PKG	32000	400	13.06	2019/11/21	A523-19110244	776780007	France	Compresse Gaze
96	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	01174T2	PKG	80000	1000	32.66	2020/1/6	A523-20010019	778210005	France	Compresse Gaze
97	Laboratoire Tetra Medical	77969/77970	Gauze Swabs with X-ray Thread	01174T2	PKG	23520	294	9.602	2020/2/22	A523-20020064	779700008	France	Compresse Gaze
98	Laboratoire Tetra Medical	77875/78183/78317	Gauze Swabs with X-ray Thread	01174T2	PKG	52000	650	21.23	2020/4/21	A523-20040496	783170002	France	Compresse Gaze
99	Laboratoire Tetra Medical	78195/78196/78305/78309/78310/	Gauze Swabs with X-ray Thread	01174T2	PKG	18000	225	7.349	2020/5/16	A523-20050251	783170002	France	Compresse Gaze
100	Laboratoire Tetra Medical	78647	Gauze Swabs with X-ray Thread	01174T2	PKG	70560	882	28.81	2020/5/23	A523-20050394	786470003	France	Compresse Gaze
101	Laboratoire Tetra Medical	77176	Gauze Swabs with X-ray Thread	01171T2	PKG	10080	84	2.743	2019/7/19	A523-19070248	771760007	France	Compresse Gaze
102	Laboratoire Tetra Medical	77493/77354	Gauze Swabs with X-ray Thread	01171T2	PKG	10080	84	2.743	2019/9/8	A523-19090045	774930005	France	Compresse Gaze
103	Laboratoire Tetra Medical	77361	Gauze Swabs with X-ray Thread	01171T2	PKG	5040	42	1.372	2019/8/15	A523-19080128	773610007	France	Compresse Gaze
104	Laboratoire Tetra Medical	77603	Gauze Swabs with X-ray Thread	01171T2	PKG	10080	84	2.743	2019/10/14	A523-19100024	776030013	France	Compresse Gaze
105	Laboratoire Tetra Medical	77680	Gauze Swabs with X-ray Thread	01171T2	PKG	10080	84	2.743	2019/10/24	A523-19100242	776800011	France	Compresse Gaze
106	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	01171T2	PKG	10080	84	2.743	2020/1/6	A523-20010019	778210003	France	Compresse Gaze
107	Laboratoire Tetra Medical	78724/78319/77877	Gauze Swabs with X-ray Thread	01171T2	PKG	20160	168	5.487	2020/6/23	A523-20060281	787240009	France	Compresse Gaze
108	Laboratoire Tetra Medical	78952	Gauze Swabs with X-ray Thread	01171T2	PKG	10080	84	2.743	2020/7/17	A523-20070203	789520005	France	Compresse Gaze

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
109	Laboratoire Tetra Medical	77179	Gauze Swabs with X-ray Thread	20602T2	PKG	6000	240	5.045	2019/7/31	A523-19070431	771790004	France	Compresse Gaze
110	Laboratoire Tetra Medical	77493/77354	Gauze Swabs with X-ray Thread	20602T2	PKG	7500	300	6.306	2019/9/8	A523-19090045	774930006	France	Compresse Gaze
111	Laboratoire Tetra Medical	77602	Gauze Swabs with X-ray Thread	20602T2	PKG	6000	240	5.045	2019/10/5	A523-19090342	776020012	France	Compresse Gaze
112	Laboratoire Tetra Medical	77683	Gauze Swabs with X-ray Thread	20602T2	PKG	1500	60	1.261	2019/11/14	A523-19110012	776830010	France	Compresse Gaze
113	Laboratoire Tetra Medical	77821	Gauze Swabs with X-ray Thread	20602T2	PKG	4500	180	3.784	2020/1/6	A523-20010019	778210007	France	Compresse Gaze
114	Laboratoire Tetra Medical	78726	Gauze Swabs with X-ray Thread	20602T2	PKG	3000	120	2.523	2020/6/16	A523-20060284	787230006	France	Compresse Gaze
115	Laboratoire Tetra Medical	78628	Gauze Swabs with X-ray Thread	20602T2	PKG	9000	360	7.568	2020/5/10	A523-20050159	786280007	France	Compresse Gaze
116	Laboratoire Tetra Medical	77352	Gauze Swabs with X-ray Thread	01029T2	PKG	10000	100	1.401	2019/8/22	A523-19080203	773520006	France	Compresse Gaze
117	Laboratoire Tetra Medical	77685	Gauze Swabs with X-ray Thread	01166T2	PKG	600	2	0.06	2019/12/13	A523-19120172	776850009	France	Compresse Gaze
118	Medline Int'l France SAS	PAN0214193	Gauze Swabs with X-ray Thread	AL017XR	CS	120	72	6.72	2019-05-08	B523-19040287	ALD19W17	France	Gauze Sponges
119	Medline Int'l France SAS	PAN0214194/ PAN0228193	Gauze Swabs with X-ray Thread	AL002XR	CS	15	120	0.795	2019-05-13	B523-19040288	ALD19W17	France	Gauze Sponges
120	Medline Int'l France SAS	PAN0214194/ PAN0228193	Gauze Swabs with X-ray Thread	AL007XR	CS	2	240	0.052	2019-05-13	B523-19040288	ALD19W17	France	Gauze Sponges
121	Medline Int'l France SAS	PAN0228191/ PAN0314195	Gauze Swabs with X-ray Thread	AL007XR	CS	30	240	0.78	2019-05-17	B523-19040398	ALD19W18	France	Gauze Sponges
122	Medline Int'l France SAS	PAN0228191/ PAN0314195	Gauze Swabs with X-ray Thread	AL021XR	CS	30	72	1.11	2019-05-17	B523-19040398	ALD19W18	France	Gauze Sponges
123	Medline Int'l France SAS	PAN0314195/ PAN0228191	Gauze Swabs with X-ray Thread	AL002XR	CS	80	120	4.24	2019-05-28	B523-19050128	ALD19W18	France	Gauze Sponges
124	Medline Int'l France SAS	PAN0314196/ PAN0327191	Gauze Swabs with X-ray Thread	AL002XR	CS	16	120	0.848	2019-06-04	B523-19050166	ALD19W20	France	Gauze Sponges

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
125	Medline Int'l France SAS	PAN0417195/ PAN0425197	Gauze Swabs with X-ray Thread	AL005XR	CS	17	120	0.952	2019-06-25	B523-19060049	ALD19W23	France	Gauze Sponges
126	Medline Int'l France SAS	PAN0417195/ PAN0425197	Gauze Swabs with X-ray Thread	AL007XR	CS	6	240	0.156	2019-06-25	B523-19060049	ALD19W23	France	Gauze Sponges
127	Medline Int'l France SAS	PAN0417195/ PAN0425197	Gauze Swabs with X-ray Thread	AL021XR	CS	8	72	0.296	2019-06-25	B523-19060049	ALD19W23	France	Gauze Sponges
128	Medline Int'l France SAS	PAN0417195/ PAN0425197	Gauze Swabs with X-ray Thread	AL002XR	CS	8	120	0.424	2019-06-25	B523-19060049	ALD19W23	France	Gauze Sponges
129	Medline Int'l France SAS	PAN0515191/ PAN0515192	Gauze Swabs with X-ray Thread	AL002XR	CS	80	120	4.24	2019-07-20	B523-19070087	ALD19W27	France	Gauze Sponges
130	Medline Int'l France SAS	PAN0515197/ PAN0417195	Gauze Swabs with X-ray Thread	AL002XR	CS	38	120	0	2019-07-20	B523-19070086	ALD19W27	France	Gauze Sponges
131	Medline Int'l France SAS	PAN0515197/ PAN0417195	Gauze Swabs with X-ray Thread	AL007XR	CS	6	240	0	2019-07-20	B523-19070086	ALD19W27	France	Gauze Sponges
132	Medline Int'l France SAS	PAN0515192/ ..PAN052719L	Gauze Swabs with X-ray Thread	AL008XR	CS	36	300	0.648	2019-07-25	B523-19070088	ALD19W28	France	Gauze Sponges
133	Medline Int'l France SAS	PAN0515192/ ..PAN052719L	Gauze Swabs with X-ray Thread	AL006XR	CS	25	120	0.875	2019-07-25	B523-19070088	ALD19W28	France	Gauze Sponges
134	Medline Int'l France SAS	PAN0515192/ ..PAN052719L	Gauze Swabs with X-ray Thread	AL017XR	CS	24	72	1.344	2019-07-25	B523-19070088	ALD19W28	France	Gauze Sponges
135	Medline Int'l France SAS	PAN0612192/ PAN0515197	Gauze Swabs with X-ray Thread	AL002XR	CS	30	120	1.59	2019-08-15	B523-19080124	ALD19W30	France	Gauze Sponges
136	Medline Int'l France SAS	PAN0612192/ PAN0515197	Gauze Swabs with X-ray Thread	AL008XR	CS	4	300	0.072	2019-08-15	B523-19080124	ALD19W30	France	Gauze Sponges
137	Medline Int'l France SAS	PAN0612192/ PAN0515197	Gauze Swabs with X-ray Thread	AL007XR	CS	8	240	0.208	2019-08-15	B523-19080124	ALD19W30	France	Gauze Sponges
138	Medline Int'l France SAS	PAN0612191/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL007XR	CS	30	240	0.78	2019-08-17	B523-19080015	ALD19W30	France	Gauze Sponges
139	Medline Int'l France SAS	PAN0612191/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL005XR	CS	56	120	3.136	2019-08-17	B523-19080015	ALD19W30	France	Gauze Sponges
140	Medline Int'l France SAS	PAN0612191/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL002XR	CS	80	120	4.24	2019-08-17	B523-19080015	ALD19W30	France	Gauze Sponges
141	Medline Int'l France SAS	PAN0612191/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL017XR	CS	30	72	1.68	2019-08-17	B523-19080015	ALD19W30	France	Gauze Sponges
142	Medline Int'l France SAS	PAN0612191/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	8	3960	0	2019-08-17	B523-19080015	ALD19W31	France	Gauze Sponges
143	Medline Int'l France SAS	PAN0627192/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL002XR	CS	80	120	4.24	2019-08-29	B523-19080225	ALD19W33	France	Gauze Sponges

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
144	Medline Int'l France SAS	PAN0627192/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL005XR	CS	112	120	6.272	2019-08-29	B523-19080225	ALD19W33	France	Gauze Sponges
145	Medline Int'l France SAS	PAN0627192/ ..PAN071519A	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	21	3960	0	2019-08-29	B523-19080225	ALD19W31	France	Gauze Sponges
146	Medline Int'l France SAS	PAN0627192	Gauze Swabs with X-ray Thread	AL026XR	CS	10	240	0.307	2019-09-07	B523-19080296	ALD19W33	France	Gauze Sponges
147	Medline Int'l France SAS	PAN062719P/ ..PAN0612192	Gauze Swabs with X-ray Thread	AL007XR	CS	6	240	0.156	2019-09-07	B523-19080295	ALD19W34	France	Gauze Sponges
148	Medline Int'l France SAS	PAN062719P/ ..PAN0612192	Gauze Swabs with X-ray Thread	AL008XR	CS	5	300	0.09	2019-09-07	B523-19080295	ALD19W34	France	Gauze Sponges
149	Medline Int'l France SAS	PAN062719P/ ..PAN0612192	Gauze Swabs with X-ray Thread	AL005XR	CS	15	120	0.84	2019-09-07	B523-19080295	ALD19W34	France	Gauze Sponges
150	Medline Int'l France SAS	PAN062719P/ ..PAN0612192	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	4	3960	3.6	2019-09-07	B523-19080295	ALD19W31	France	Gauze Sponges
151	Medline Int'l France SAS	PAN0716192/ PAN0627192	Gauze Swabs with X-ray Thread	AL017XR	CS	120	72	6.72	2019-09-24	B523-19090237	ALD19W33	France	Gauze Sponges
152	Medline Int'l France SAS	PAN0716191/ PAN062719P	Gauze Swabs with X-ray Thread	AL005XR	CS	15	120	0.84	2019-09-24	B523-19090045	ALD19W36	France	Gauze Sponges
153	Medline Int'l France SAS	PAN0716191/ PAN062719P	Gauze Swabs with X-ray Thread	AL021XR	CS	3	72	0.111	2019-09-24	B523-19090045	ALD19W36	France	Gauze Sponges
154	Medline Int'l France SAS	PAN0716191/ PAN062719P	Gauze Swabs with X-ray Thread	AL002XR	CS	26	120	1.378	2019-09-24	B523-19090045	ALD19W36	France	Gauze Sponges
155	Medline Int'l France SAS	PAN0716191/ PAN062719P	Gauze Swabs with X-ray Thread	AL007XR	CS	4	240	0.104	2019-09-24	B523-19090045	ALD19W36	France	Gauze Sponges
156	Medline Int'l France SAS	PAN0716192/ PAN0725192	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	21	3960	14.71	2019-09-30	B523-19090238	ALD19W37	France	Gauze Sponges
157	Medline Int'l France SAS	PAN0716192/ PAN0725192	Gauze Swabs with X-ray Thread	AL005XR	CS	56	120	3.136	2019-09-30	B523-19090238	ALD19W37	France	Gauze Sponges
158	Medline Int'l France SAS	PAN0716192/ PAN0725192	Gauze Swabs with X-ray Thread	AL002XR	CS	40	120	2.12	2019-09-30	B523-19090238	ALD19W39	France	Gauze Sponges
159	Medline Int'l France SAS	PAN0716192/ PAN0725192	Gauze Swabs with X-ray Thread	AL017XR	CS	40	72	2.24	2019-09-30	B523-19090238	ALD19W39	France	Gauze Sponges
160	Medline Int'l France SAS	PAN0725192	Gauze Swabs with X-ray Thread	AL007XR	CS	25	240	0.65	2019-10-11	B523-19090235	ALD19W39	France	Gauze Sponges
161	Medline Int'l France SAS	PAN0814192	Gauze Swabs with X-ray Thread	AL007XR	CS	4	240	0.104	2019-10-17	B523-19090362	ALD19W40	France	Gauze Sponges
162	Medline Int'l France SAS	PAN0814192	Gauze Swabs with X-ray Thread	AL008XR	CS	5	300	0.09	2019-10-17	B523-19090362	ALD19W40	France	Gauze Sponges

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
163	Medline Int'l France SAS	PAN0814192	Gauze Swabs with X-ray Thread	AL002XR	CS	24	120	1.272	2019-10-17	B523-19090362	ALD19W40	France	Gauze Sponges
164	Medline Int'l France SAS	PAN0814192	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	13	3960	11.7	2019-10-17	B523-19090362	ALD19W40	France	Gauze Sponges
165	Medline Int'l France SAS	PAN0814192	Gauze Swabs with X-ray Thread	AL021XR	CS	4	72	0.148	2019-10-17	B523-19090362	ALD19W40	France	Gauze Sponges
166	Medline Int'l France SAS	PAN0814191/ ..PAN0919191	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	9	3960	8.1	2019-10-29	B523-19100087	ALD19W41	France	Gauze Sponges
167	Medline Int'l France SAS	PAN0814191/ ..PAN0919191	Gauze Swabs with X-ray Thread	AL005XR	CS	56	120	3.136	2019-10-29	B523-19100087	ALD19W41	France	Gauze Sponges
168	Medline Int'l France SAS	PAN0919198/ PAN0814192	Gauze Swabs with X-ray Thread	AL008XR	CS	2	300	0.036	2019-11-17	B523-19100426	ALD19W44	France	Gauze Sponges
169	Medline Int'l France SAS	PAN0919198/ PAN0814192	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	8	3960	7.2	2019-11-17	B523-19100426	ALD19W44	France	Gauze Sponges
170	Medline Int'l France SAS	PAN0919198/ PAN0814192	Gauze Swabs with X-ray Thread	AL007XR	CS	4	240	0.104	2019-11-17	B523-19100426	ALD19W44	France	Gauze Sponges
171	Medline Int'l France SAS	PAN0919198/ PAN0814192	Gauze Swabs with X-ray Thread	AL005XR	CS	18	120	1.008	2019-11-17	B523-19100426	ALD19W44	France	Gauze Sponges
172	Medline Int'l France SAS	PAN0919198/ PAN0814192	Gauze Swabs with X-ray Thread	AL021XR	CS	2	72	0.074	2019-11-17	B523-19100426	ALD19W44	France	Gauze Sponges
173	Medline Int'l France SAS	PAN0919199/ PAN0919198	Gauze Swabs with X-ray Thread	AL021XR	CS	3	72	0.111	2019-11-22	B523-19110079	ALD19W45	France	Gauze Sponges
174	Medline Int'l France SAS	PAN0919199/ PAN0919198	Gauze Swabs with X-ray Thread	AL005XR	CS	18	120	1.008	2019-11-22	B523-19110079	ALD19W45	France	Gauze Sponges
175	Medline Int'l France SAS	PAN0919199/ PAN0919198	Gauze Swabs with X-ray Thread	AL002XR	CS	24	120	1.272	2019-11-22	B523-19110079	ALD19W45	France	Gauze Sponges
176	Medline Int'l France SAS	PAN0919199/ PAN0919198	Gauze Swabs with X-ray Thread	AL007XR	CS	7	240	0.182	2019-11-22	B523-19110079	ALD19W45	France	Gauze Sponges
177	Medline Int'l France SAS	PAN0919199/ PAN0919198	Gauze Swabs with X-ray Thread	AL008XR	CS	4	300	0.072	2019-11-22	B523-19110079	ALD19W45	France	Gauze Sponges
178	Medline Int'l France SAS	PAN0919192	Gauze Swabs with X-ray Thread	AL002XR	CS	40	120	2.12	2019-11-24	B523-19110078	ALD19W45	France	Gauze Sponges
179	Medline Int'l France SAS	PAN0919191	Gauze Swabs with X-ray Thread	AL017XR	CS	40	72	2.24	2019-11-24	B523-19100311	ALD19W43	France	Gauze Sponges
180	Medline Int'l France SAS	PAN0919191	Gauze Swabs with X-ray Thread	AL005XR	CS	28	120	1.568	2019-11-24	B523-19100311	ALD19W43	France	Gauze Sponges
181	Medline Int'l France SAS	PAN0919192	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	19	3960	17.1	2019-11-24	B523-19110078	ALD19W45	France	Gauze Sponges

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
182	Medline Int'l France SAS	PAN0919192	Gauze Swabs with X-ray Thread	AL006XR	CS	25	120	0.875	2019-11-24	B523-19110078	ALD19W45	France	Gauze Sponges
183	Medline Int'l France SAS	PAN0919192	Gauze Swabs with X-ray Thread	AL017XR	CS	40	72	2.24	2019-11-24	B523-19110078	ALD19W45	France	Gauze Sponges
184	Medline Int'l France SAS	PAN0919193	Gauze Swabs with X-ray Thread	AL008XR	CS	25	300	0.45	2019-11-27	B523-19110167	ALD19W46	France	Gauze Sponges
185	Medline Int'l France SAS	PAN0919193	Gauze Swabs with X-ray Thread	AL017XR	CS	40	72	2.24	2019-11-27	B523-19110167	ALD19W46	France	Gauze Sponges
186	Medline Int'l France SAS	PAN0919193	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	11	3960	9.9	2019-11-27	B523-19110167	ALD19W46	France	Gauze Sponges
187	Medline Int'l France SAS	PAN0919193	Gauze Swabs with X-ray Thread	AL005XR	CS	56	120	3.136	2019-11-27	B523-19110167	ALD19W46	France	Gauze Sponges
188	Medline Int'l France SAS	PAN0919193	Gauze Swabs with X-ray Thread	AL002XR	CS	40	120	2.12	2019-11-27	B523-19110167	ALD19W46	France	Gauze Sponges
189	Medline Int'l France SAS	PAN0919199	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	5	3960	4.5	2019-12-02	B523-19120019	ALD19W45	France	Gauze Sponges
190	Medline Int'l France SAS	PAN0919194/ ..PAN0919193	Gauze Swabs with X-ray Thread	AL017XR	CS	40	72	2.24	2019-12-07	B523-19110263	ALD19W47	France	Gauze Sponges
191	Medline Int'l France SAS	PAN0919194	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	4	3960	3.6	2019-12-10	B523-19120081	ALD19W47	France	Gauze Sponges
192	Medline Int'l France SAS	PAN091919A	Gauze Swabs with X-ray Thread	AL021XR	CS	1	72	0.037	2019-12-13	B523-19110381	ALD19W48	France	Gauze Sponges
193	Medline Int'l France SAS	PAN091919A	Gauze Swabs with X-ray Thread	AL008XR	CS	4	300	0.072	2019-12-13	B523-19110381	ALD19W48	France	Gauze Sponges
194	Medline Int'l France SAS	PAN091919A	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	7	3960	6.3	2019-12-13	B523-19110381	ALD19W48	France	Gauze Sponges
195	Medline Int'l France SAS	PAN091919A	Gauze Swabs with X-ray Thread	AL007XR	CS	6	240	0.156	2019-12-13	B523-19110381	ALD19W48	France	Gauze Sponges
196	Medline Int'l France SAS	PAN091919A	Gauze Swabs with X-ray Thread	AL002XR	CS	24	120	1.272	2019-12-13	B523-19110381	ALD19W48	France	Gauze Sponges
197	Medline Int'l France SAS	PAN091919A	Gauze Swabs with X-ray Thread	AL005XR	CS	18	120	1.008	2019-12-13	B523-19110381	ALD19W48	France	Gauze Sponges
198	Medline Int'l France SAS	PAN0919195/ PAN0919194	Gauze Swabs with X-ray Thread	AL017XR	CS	80	72	4.48	2019-12-17	B523-19110382	ALD19W48	France	Gauze Sponges
199	Medline Int'l France SAS	PAN0919195/ PAN0919194	Gauze Swabs with X-ray Thread	AL007XR	CS	25	240	0.65	2019-12-17	B523-19110382	ALD19W48	France	Gauze Sponges
200	Medline Int'l France SAS	PAN0919195/ PAN0919194	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	9	3960	8.1	2019-12-17	B523-19110382	ALD19W48	France	Gauze Sponges

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
201	Medline Int'l France SAS	PAN0919195/ PAN0919194	Gauze Swabs with X-ray Thread	AL005XR	CS	56	120	3.136	2019-12-17	B523-19110382	ALD19W48	France	Gauze Sponges
202	Medline Int'l France SAS	PAN1017192/ PAN091919A	Gauze Swabs with X-ray Thread	AL005XR	CS	18	120	1.008	2020-01-13	B523-19120140	ALD19W50	France	Gauze Sponges
203	Medline Int'l France SAS	PAN1017192/ PAN091919A	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	1	3960	0.9	2020-01-13	B523-19120140	ALD19W50	France	Gauze Sponges
204	Medline Int'l France SAS	PAN1022191/ ..PAN0919196	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	12	3960	10.8	2020-01-15	B523-19120213	ALD19W52	France	Gauze Sponges
205	Medline Int'l France SAS	PAN1022191/ ..PAN0919196	Gauze Swabs with X-ray Thread	AL002XR	CS	40	120	2.12	2020-01-15	B523-19120213	ALD19W52	France	Gauze Sponges
206	Medline Int'l France SAS	PAN1030191/ ..PAN1017191	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	4	3960	3.6	2020-01-19	B523-19120363	ALD19W53	France	Gauze Sponges
207	Medline Int'l France SAS	PAN1030191/ ..PAN1017191	Gauze Swabs with X-ray Thread	AL017XR	CS	40	72	2.24	2020-01-19	B523-19120363	ALD19W53	France	Gauze Sponges
208	Medline Int'l France SAS	PAN1112193/ ..PAN122019X	Gauze Swabs with X-ray Thread	AL017XR	CS	80	72	4.492	2020-03-05	B523-20020067	ALD20W03	France	Gauze Sponges
209	Medline Int'l France SAS	PAN1112192/ ..PAN1217195	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	2	3960	1.8	2020-03-17	B523-20030115	ALD20W06	France	Gauze Sponges
210	Medline Int'l France SAS	PAN1112192/ ..PAN1217195	Gauze Swabs with X-ray Thread	AL005XR	CS	18	120	1.008	2020-03-17	B523-20030115	ALD20W01	France	Gauze Sponges
211	Medline Int'l France SAS	PAN1112192/ ..PAN1217195	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	5	3960	4.5	2020-03-17	B523-20030115	ALD20W07	France	Gauze Sponges
212	Medline Int'l France SAS	PAN1112192/ ..PAN1217195	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	7	3960	6.3	2020-03-27	B523-20030116	ALD20W07	France	Gauze Sponges
213	Medline Int'l France SAS	PAN011620F/ ..PAN0305203	Gauze Swabs with X-ray Thread	AL002XR	CS	25	120	1.336	2020-05-07	B523-20040349	ALD20W17	France	Gauze Sponges
214	Medline Int'l France SAS	PAN0116203/ PAN020620T	Gauze Swabs with X-ray Thread	AL002XR	CS	40	120	2.137	2020-05-19	B523-20050121	ALD20W18	France	Gauze Sponges
215	Medline Int'l France SAS	PAN0116203/ PAN020620T	Gauze Swabs with X-ray Thread	AL003XR- BULK	CS	1	3960	0.996	2020-05-19	B523-20050121	ALD20W18	France	Gauze Sponges
216	Medline Int'l France SAS	PAN020620U/ ..PAN1217191	Gauze Swabs with X-ray Thread	AL001XR- BULK	CS	3	1980	2.7	2020-05-19	B523-20040353	ALD20W19	France	Gauze Sponges
217	Medline Int'l France SAS	PAN011620F/ ..PAN021120A	Gauze Swabs with X-ray Thread	AL021XR	CS	25	72	0.915	2020-05-22	B523-20050061	ALD20W19	France	Gauze Sponges
218	Medline Int'l France SAS	PAN011620F/ ..PAN021120A	Gauze Swabs with X-ray Thread	AL002XR	CS	24	120	1.282	2020-05-22	B523-20050061	ALD20W19	France	Gauze Sponges
219	Medline Int'l France SAS	PAN011620F/ ..PAN021120A	Gauze Swabs with X-ray Thread	AL001XR- BULK	CS	11	1980	9.9	2020-05-22	B523-20050061	ALD20W17	France	Gauze Sponges

Item	Customer	Purchase No.	Product Category	Product Code.	Unit	Qty Pcs	Qty Box	Volume	On Board Date	Invoice No.	Batch No.	End Market	Trade Name
220	Medline Int'l France SAS	PAN011620F/ ..PAN020620S	Gauze Swabs with X-ray Thread	AL002XR	CS	25	120	1.336	2020-05-22	B523-20050062	ALD20W19	France	Gauze Sponges
221	Medline Int'l France SAS	PAN0305204/ ..PAN041020G	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	5	3960	4.5	2020-06-29	B523-20060135	ALD20W19	France	Gauze Sponges
222	Medline Int'l France SAS	PAN0305204/ ..PAN041020G	Gauze Swabs with X-ray Thread	AL026XR	CS	10	240	0.307	2020-06-29	B523-20060135	ALD20W24	France	Gauze Sponges
223	Medline Int'l France SAS	PAN041020H	Gauze Swabs with X-ray Thread	AL003XR- BULK	CS	13	3960	12.95	2020-06-30	B523-20060137	ALD20W24	France	Gauze Sponges
224	Medline Int'l France SAS	PAN050820P/ PAN052120Q	Gauze Swabs with X-ray Thread	AL008XR	CS	25	300	0.447	2020-07-24	B523-20070093	ALD20W27	France	Gauze Sponges
225	Medline Int'l France SAS	PAN050820P/ PAN052120Q	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	5	3960	4.5	2020-07-24	B523-20070093	ALD20W27	France	Gauze Sponges
226	Medline Int'l France SAS	PAN050820R/ ..PAN052920A	Gauze Swabs with X-ray Thread	AL005XR	CS	25	120	1.389	2020-07-24	B523-20070095	ALD20W28	France	Gauze Sponges
227	Medline Int'l France SAS	PAN050820R/ ..PAN052920A	Gauze Swabs with X-ray Thread	AL002XR	CS	25	120	1.336	2020-07-24	B523-20070095	ALD20W28	France	Gauze Sponges
228	Medline Int'l France SAS	PAN050820Q/ ..PAN041020P	Gauze Swabs with X-ray Thread	AL002XR	CS	80	120	4.274	2020-07-28	B523-20070094	ALD20W28	France	Gauze Sponges
229	Medline Int'l France SAS	PAN050820Q/ ..PAN041020P	Gauze Swabs with X-ray Thread	AL004XR A-BULK	CS	20	3960	18	2020-07-28	B523-20070094	ALD20W27	France	Gauze Sponges

Remark:

1. Due to we expand the scope of internal retained sample batches confirm, in the meanwhile we were keeping communication with our customer and collecting information from market, so we updated the form information as above to secure fields influence.
2. The information might be change due to the reason above.