

TABLE DE CODAGE hGH - IGF-I - novembre 2016

(liste non exhaustive)

PRINCIPAUX AUTOMATES

BECKMAN COULTER	Access / Access 2	ULA
	Lxi 725	DCM
	UniCel DxI 800	UCD
	UniCel Dx C 600i	DCP
DIASORIN S.A.	Liaison XL	UKV
	Liaison	UKW
ORTHO Clin. Diag.	Vitros ECI	U4V
	Vitros 3600	U4W
	Vitros 5600	FKI
PERKIN ELMER Inst	DELFI A Xpress	U4A
	AutoDELFIA	U4D
ROCHE Diagnostics	Elecsys 2010	UWG
	Modular E (170)/EE	UWH
	Cobas 4000 e411	UWL
	Cobas 6000 e601	UWR
	Cobas 8000 e602	UWT
	Cobas 8000 e801	UWS
SIEMENS MED. SOL. DIAG.	ADVIA Centaur / ADVIA Centaur XP / ADVIA Centaur XPT	U4S
	ADVIA Centaur CP	U4E
	Dimension ARx/XL/RxL	DFG
	Dimension Xpand	DFI
	Dimension RxL HM (Heterogeneous Module)	DFH
	Dimension Vista	DFJ
	Dimension ExL	DFK
	Immolute / Immolute 1000	U4L
	Immolute 2000 / Immolute 2000XPi	U4R
	Immolute 2500	U4P
TOSOH BIOSCIENCE	AIA-600	UEB
	AIA-21 / AIA-600-II	UEC
	AIA 360	UEN
	AIA 1800	UEP
	AIA 2000	UER
	AIA-900	UET
Tout lecteur de plaque		ZPX
Tout compteur gamma pour dosage radioimmunologique		5CX
Autres		XXX

TABLE DE CODAGE hGH - IGF-I - novembre 2016 - (Liste non exhaustive)
CODE TECHNIQUES - REACTIFS

DISTRIBUTEUR	NOM DU COFFRET	REFERENCE	CODE
hGH			
BECKMAN COULTER	Access Ultrasensitive hGH	33580	QE
	IRMA GH (Immunotech)	IM 1397	BO
EUROIMMUN - BIO-ADVANCE	DRG ELISA hGH	EIA 1787	6B
Cisbio Bioassays	hGH-RIACT	hGH-RIACT	BN
DIASORIN S.A.	Liaison hGH	310340	S8
DIASOURCE	hGH IRMA	KIP 1081	AM
IDS	hGH iSYS	IS-3700	SC
PERKIN ELMER	DELFI/AutoDELFIA hGH	1244 - 041 / B 041-101	KC
ROCHE Diag.	ELECSYS hGH	5390125190	RD
SIEMENS MED.SOL. DIAGN.	Immulin / Im. 2000-2000 Xpi hGH	6602230 / 6606147	SA
TOSOH Biosciences	Stat AIA Pack hGH	025266	DL
IGF-I			
BECKMAN COULTER	IRMA IGF-I	A15729	AZ
EUROIMMUN - BIO-ADVANCE	ELISA IGF-1 600	EIA 4140	MB
Cisbio Bioassays	Insulin-like Growth Factor I (Somatomedin C)	IGF1-RIACT	BN
DIASORIN S.A.	LIAISON® IGF-I	313231	S8
	IGF-1 RIA-CT (Mediagnost)	IGF-R22	B9
DIASOURCE	RIA-CT-IGF-1	KIP-1588	AB
IDS	IDS-iSYS Insulin like Growth Factor-I	IS-3900	SC
SIEMENS MED.SOL. DIAGN.	IMMULITE/ IM. 2000/ 2000 XPi IGF-1 ancienne trousse	LKGF1/LKGF5/ LKGF2	SA
	IMMULITE/ IM. 2000/ 2000 XPi IGF-1 nouvelle trousse	LKIGF1/ L2KIGF2	SB