

CONTRÔLE NATIONAL DE QUALITE TABLES DE CODAGE BIOCHIMIE 15B10-1

(mars 2015)

Liste non exhaustive

Si vous ne trouvez pas le réactif et/ou l'appareil utilisé, indiquez en clair, sur le formulaire de réponse, les informations sur le réactif (distributeur, nom du coffret, référence catalogue) et/ou l'automate (distributeur, fabricant, modèle...). Si possible, nous adresser par courrier ou par mail (cnq.labm@ansm.sante.fr), la fiche technique (notice) du réactif.

TROPONINE I

Code	Distributeur ou fabricant	Nom du coffret ou appareil utilisé
RJ	ABBOTT	, ARCHITECT i-systems STAT Troponin-I - REF. 2K41
RJ	ABBOTT	, ARCHITECT i-systems STAT hs-Troponin-I - REF. 3P25
DJ	ABBOTT	, AxSYM Troponin-I ADV - REF. 2J44
(*)	ABBOTT	, i-STAT cTnI - REF. 03P90-25
QE	BECKMAN COULTER	, Access-TnI - REF. A78803
QE	BECKMAN COULTER	, Access AccuTnI+3 w/Access 2, DxI platforms - REF. A98264
DB	BIOMERIEUX	, VIDAS Troponin I Ultra - REF. 30448
ZI	BIOSITE (ALERE)	, Triage Cardiac Panel (TnI) - REF. 97400, 97500, 97600
S8	DIASORIN	, LIAISON Troponin I - REF. 315101
(*)	MITSUBISHI	, PATHFAST cTnI - REF. 1110-2000
P5	ORTHO-CD	, VITROS Troponin I ES - REF. 680 2301
KR	RADIOMETER	, AQT90 FLEX TnI - REF. 942-903
(*)	RESPONSE Biomedical	, RAMP Troponin I
RD	ROCHE	, Elecsys Troponin I - REF. 05094810
RD	ROCHE	, Elecsys Troponin I STAT - REF. 05094798
SI	SIEMENS	, ADVIA Centaur TnI Ultra - REF. 2789602, 2790309
NA	SIEMENS	, Dimension RxL/Xpand CTNI - Troponin-Ic - REF. RF421C
NA	SIEMENS	, Dimension RxL/Xpand LTNI - Troponin-Ic - REF. RF521
SQ	SIEMENS	, Dimension EXL TNI - LOCI Troponin-Ic - REF. RF621
SQ	SIEMENS	, Dimension Vista CTNI - LOCI Troponin-Ic - REF. K6421
SA	SIEMENS	, IMMULITE Troponin I - REF. 6603512
SA	SIEMENS	, IMMULITE Turbo Troponin I - REF. 6604454
SA	SIEMENS	, IMMULITE 2000 Troponin I - REF. 6607445
E2	SIEMENS	, Stratus CS Acute Care Troponin-I TestPak - REF. CCTNI
DL	TOSOH	, AIA series ST AIA-PACK cTnI 2nd-Gen (cTnI_2) - REF. 0025205
DL	TOSOH	, AIA series ST AIA-PACK cTnI 3rd-Gen (cTnI_3) - REF. 0025215

TROPONINE T

Code	Distributeur ou fabricant	Nom du coffret ou appareil utilisé
KR	RADIOMETER	, AQT90 FLEX TnT - REF. 942-940
RD	ROCHE	, Elecsys Troponin T hs - REF. 05092744
RD	ROCHE	, Elecsys Troponin T hs STAT - REF. 05092728
(*)	ROCHE	, Cardiac Reader Cardiac Troponin T Quantitative - REF. 11894307
(*)	ROCHE	, cobas h-232 Cardiac Troponin T Quantitative - REF. 04877772

HEMOGLOBINE A1c (HbA1c)

Code	Distributeur ou fabricant	Nom du coffret ou appareil utilisé
9J	ABBOTT	, ARCHITECT c-systems Hemoglobin A1c - REF. 4P52
RJ	ABBOTT	, ARCHITECT i-systems Hemoglobin A1c - REF. 4P72
H6	BECKMAN COULTER	, AU systems HbA1c - REF. OSR6192
H6	BECKMAN COULTER	, AU systems HbA1c - REF. B00389
HK	BECKMAN COULTER	, DxC 600/800 HbA1c- - REF. A87905
HK	BECKMAN COULTER	, DxC 600/800 HbA1c3 - REF. B36415
HS	BIOLABO	, HbA1c - REF. 22010, 22011
20	BIO-RAD	, VARIANT II Hemoglobin A1c
2M	BIO-RAD	, VARIANT II TURBO Hemoglobin A1c
2P	BIO-RAD	, D-10 HbA1c
HT	DIASYS	, One HbA1c FS - REF. 1 3329 ..
(*)	DIASYS	, HbA1c net FS - REF. 1 3348 ..
2G	ELITECH (ARKRAY)	, ADAMS A1c HA-8180V
HW	HORIBA ABX	, HbA1c WB - REF. A11A01702
2F	MENARINI (ARKRAY)	, ADAMS A1c HA-8140
2C	MENARINI (ARKRAY)	, ADAMS A1c HA-8160
(*)	MENARINI (TRINITY Biotech)	, Premier Hb9210 HbA1c
HF	ORTHO-CD	, VITROS d%A1c - REF. 680 2314
AN	PROGEN (AXIS-SHIELD)	, NycoCard HbA1c
HH	RANDOX	, HbA1c - REF. 3443
HO	ROCHE	, Hitachi/Modular P HbA1c Gen.3 - REF. 05336104, 05335787
HZ	ROCHE	, Integra 400/800 HbA1c Gen.2 (A1C-2) - REF. 04528123
HR	ROCHE	, cobas c-311, c-501/2 & c-701/2 HbA1c Gen.3 (A1C-3) - REF. 05336163
(*)	ROCHE	, cobas b-101 POC system - REF. 06378676, 06777201
1C	SEBIA	, CAPILLARYS Hb A1c w/Capillarys 2 Flex Piercing
(*)	SEBIA	, CAPI 3 Hb A1c - REF. 2515 w/Capillarys 3
1M	SEBIA	, MINICAP Hb A1c w/Minicap Flex Piercing
HE	SIEMENS	, ADVIA Chemistry systems HbA1C-3 - REF. 10379673, 10485591
HA	SIEMENS	, Dimension RxL/Xpand/EXL HB1C - REF. DF105A
HA	SIEMENS	, Dimension Vista HA1C - REF. K3105A
HB	SIEMENS	, DCA 2000+/DCA Vantage - REF. 6162000
HY	THERMO Scientific	, HbA1c w/Konelab, Indiko - REF. 981947, 981948
2T	TOSOH	, G7
2R	TOSOH	, G8
(*)	TOSOH	, GX (HLC-723GX)

APPAREILS

Code	Distributeur ou fabricant	Modèle
DZG	ABBOTT	ARCHITECT c4000
DZH	ABBOTT	ARCHITECT c8000
DZI	ABBOTT	ARCHITECT c16000
U4Y	ABBOTT	ARCHITECT i1000SR / i2000SR
UOF	ABBOTT	AxSYM/AxSYM +
WPA	ABBOTT	i-STAT analyzer
DEL	BECKMAN COULTER	AU680
DEP	BECKMAN COULTER	AU480
DER	BECKMAN COULTER	AU400
DES	BECKMAN COULTER	AU640
DET	BECKMAN COULTER	AU2700

DEW	BECKMAN COULTER	AU5400
DEZ	BECKMAN COULTER	AU5800
DCM	BECKMAN COULTER	Synchron LXi 725
DCP	BECKMAN COULTER	UniCel DxC 600/600i
DCQ	BECKMAN COULTER	UniCel DxC 800
UCD	BECKMAN COULTER	UniCel DxI 600/800
ULA	BECKMAN COULTER	Access/Access 2
UGV	BIOMERIEUX	VIDAS
UGW	BIOMERIEUX	mini VIDAS
UGT	BIOMERIEUX	VIDAS 3
9GB	BIO-RAD	D10
9GR	BIO-RAD	VARIANT & VARIANT II
9GS	BIO-RAD	VARIANT II TURBO
TRI	BIOSITE (ALERE)	Alere Triage Meters
UKV	DIASORIN	Liaison XL
UKW	DIASORIN	Liaison
9GE	ELITECH (ARKRAY)	ADAMS HA-8180
DAP	HORIBA ABX	Pentra 400
DAQ	HORIBA ABX	Pentra C200
DKM	MAXMAT S. A.	Maxmat PL
9GN	MENARINI (ARKRAY)	ADAMS HA-8140
9GP	MENARINI (ARKRAY)	ADAMS HA-8160
(*)	MENARINI (TRINITY)	Premier Hb9210
(*)	mitsubishi	Pathfast
FKG	ORTHO-CD	Vitros 5,1 FS (Fusion)
FKI	ORTHO-CD	Vitros 5600
FKJ	ORTHO-CD	Vitros 4600
U4V	ORTHO-CD	Vitros ECi/ECiQ
U4W	ORTHO-CD	Vitros 3600
9GA	PROGEN (AXIS-SHIELD)	NycoCard reader
UCR	RADIOMETER	AQT90 FLEX analyzer
(*)	RESPONSE Biomedical	RAMP Reader
(*)	ROCHE	Cardiac Reader
(*)	ROCHE	cobas b-101 POC system
(*)	ROCHE	cobas h-232
DQC	ROCHE	cobas c-111
DQH	ROCHE	Cobas Integra 400/400 +
DQL	ROCHE	Cobas Integra 800
DQN	ROCHE	cobas c-311 (cobas 4000 series)
DQP	ROCHE	cobas c-501 (cobas 6000 series)
DQR	ROCHE	cobas c-701/c-502 (cobas 8000 series)
DWE	ROCHE	Hitachi 717
DWG	ROCHE	Hitachi 911
DWI	ROCHE	Modular P/PP/DP
DWL	ROCHE	Hitachi 912
UWG	ROCHE	Elecsys 2010
UWG	ROCHE	Elecsys 2010
UWH	ROCHE	Modular E (170)/EE
UWH	ROCHE	Modular E (170)/EE
UWL	ROCHE	cobas e-411 (cobas 4000 series)
UWR	ROCHE	cobas e-601 (cobas 6000 series)
UWT	ROCHE	cobas e-602 (cobas 8000 series)

VMC	SEBIA	Capillarys/Capillarys 2/Capillarys 2 Flex Piercing
VMD	SEBIA	Minicap/Minicap Flex Piercing
(*)	SEBIA	Capillarys 3 TERA
DFH	SIEMENS	Dimension RxL/RxL Max w/HM
DFI	SIEMENS	Dimension Xpand/Xpand Plus w/HM
DFJ	SIEMENS	Dimension Vista
DFK	SIEMENS	Dimension EXL
DTL	SIEMENS	ADVIA 1200
DTM	SIEMENS	ADVIA 2400
DTN	SIEMENS	ADVIA Chemistry XPT
DTY	SIEMENS	ADVIA 1650/1800
U4E	SIEMENS	ADVIA Centaur CP
U4S	SIEMENS	ADVIA Centaur / ADVIA Centaur XP
U4L	SIEMENS	Immulite/Immulite 1000
U4R	SIEMENS	Immulite 2000/Immulite 2000 XPi
UBB	SIEMENS	DCA 2000 +
UBC	SIEMENS	DCA Vantage
UYB	SIEMENS	Stratus CS
DBL	THERMO Scientific	Konelab 30/i & PRIME 30/i
DBM	THERMO Scientific	Konelab 60/i & PRIME 60/i
DBN	THERMO Scientific	Konelab 20/i
DBQ	THERMO Scientific	Konelab 20XT/i
DBR	THERMO Scientific	Indiko/Indiko Plus
9GW	TOSOH	G7
9GY	TOSOH	G8
(*)	TOSOH	GX (HLC-723GX)
UEC	TOSOH	AIA-21/AIA-600 II
UEN	TOSOH	AIA-360
UEP	TOSOH	AIA-1800
UER	TOSOH	AIA-2000
UET	TOSOH	AIA-900

(*) Encodage en cours, le cas échéant.