

Consensus 75% anticorps anti-HLA

toutes techniques confondues

échantillon	dépistage (1)				identification (1)										(2)
15S1		I	G		B44	B45	B13	B76	B18	B27	B37	B47	B82	Cw15	25
					Cw17	Cw18	Cw5	Cw6							
		II	G		DQ5	DQ6	DQ8	DQ9							25
15S2		I	G		A66	B45	B13	B64	B65	B75	B62	B76	B39	B18	25
					B50	B56	B54	B55	B27	B35	B60	B61	B41	B42	
					B47	B48	B67	B7	B72	B71	B73	B78	B8	B81	
					B82										
	N														18
15S3		I	G		B45	B44	B76	B75							23
		II	G		DR17	DR18	DR12	DR11	DR52	DR14	DR13	DR8			25
15S4	N														19
	N														17

(1) : I : positif classe I ; II : positif classe II ; N : négatif (classe I ou classe II) ; G : IgG ; M : IgM ; NC : non consensus ; na : non applicable

(2) : nombre de laboratoires ayant rendu au moins une spécificité dans la classe

(1) : I : positif classe I ; II : positif classe II ; N : négatif (classe I ou classe II) ; G : IgG ; M : IgM ; NC : non consensus ; na : non applicable
(2) : nombre de laboratoires ayant rendu au moins une spécificité dans la classe