

**TABLE DE CODAGE GAZOMETRIE SANGUINE 2015
(mars 2015, liste non exhaustive)**

DISTRIBUTEUR	MODELE	CODE
ROCHE Diagnostics	Cobas b 121 ex Omni C	WAM
	Cobas b 221 ex Omni S	WAS
	Cobas b123	WAT
OSMETECH / MSR SEBAC DIAG.	Opti 1	WAJ
	Opti 3 - CcA	WAK
SIEMENS MED. SOL. DIAG.	Rapidpoint 400/405	WBB
	Rapidlab 348	WCM
	Rapidlab 238-248	WCO
	Rapidlab série 800/840/845	WCP
	Rapidlab série 850/855	WCQ
	Rapidlab série 860/865	WCR
	Rapidlab série 1200/1265	WCS
	Rapidpoint 340/350	WCT
	Rapidpoint 500	WCU
GAMIDA	Diametrics IRMA SL série 2000	WDA
ALERE	EPOC	WFA
INSTRUMENTATION LABORATORY	GEM 3000	WIP
	GEM 4000	WIQ
	GEM 3500	WIT
	Synthesis	WIR
NOVA Biomedical	pHOx - pHOx Ultra	WNM
	CCX1-1+CCX2-2+	WNO
	CCX3-3+	WNP
ABBOTT	I Stat	WPA
RADIOMETER	ABL 500 - 600 - 605	WRB
	ABL 555	WRD
	ABL 520 - 620 - 625	WRH
	ABL 5/ BPH 5	WRI
	ABL 70	WRL
	ABL série 700	WRM
	NPT7	WRN
	ABL série 800	WRO
	ABL 80	WRP
	ABL 90	WRS